

## Spec-10:2K

2048 x 512 imaging array | 13.5 x 13.5- $\mu$ m pixels



The Spec-10:2K from Princeton Instruments utilizes a high-performance, back-illuminated CCD with a format ideal for high-resolution spectroscopy applications. The UV/AR-coated detector offers excellent sensitivity in the ultraviolet. A "Princeton Instruments exclusive" CCD feature provides software-selectable amplifiers that permit operation in either high-capacity mode (absorbance spectroscopy) or high-sensitivity mode (Raman or fluorescence spectroscopy). This feature delivers sensitivity and dynamic range unmatched by industry-standard "1024 pixel" CCDs. Cooling the CCD to cryogenic temperatures effectively eliminates dark noise, even for long exposure times.


Features	Benefits
<b>2048 x 512 imaging array</b>	Provides highest level of resolution for demanding applications
<b>13.5 x 13.5-<math>\mu</math>m pixels</b>	Small pixel size supports high resolution
<b>Back-illuminated, UV/AR-coated CCD</b>	High quantum efficiency for low-light applications in the ultraviolet
<b>Cryogenic cooling</b>	Effective elimination of dark noise, even for long exposure times
<b>Software-selectable amplifiers</b>	Exclusive feature provides highest level of sensitivity and dynamic range for absorbance, Raman, and fluorescence applications
<b>Standard spectrometer interface</b>	Will interface with most spectrometers
<b>Dual-digitizer option</b>	Multiple-speed digitization allows complete freedom to select between "slow operation" for low noise and highest SNR (signal-to-noise ratio) or "fast operation" for rapid image acquisition
<b>"USB 2.0 interface" configuration</b>	Seamless, plug-and-play connection to PC notebooks and desktops Easy OEM integration
<b>"PCI interface" configuration</b>	Industry standard for fast, reliable data transfer
<b>WinSpec and PVCAM<sup>®</sup></b>	Offers easy-yet-sophisticated Windows <sup>®</sup> GUI controls Automates data acquisition, analysis and display
<b>Linux<sup>®</sup> drivers and SITK<sup>™</sup> plug-in for National Instruments' LabVIEW<sup>™</sup></b>	Extends system utility



## Spec-10:2K Specifications

	Spec-10:2KB	Spec-10:2KBUV	
<b>CCD Image Sensor</b>	back-illuminated	back-illuminated, enhanced UV	
	All Spec-10:2Ks		
	e2v CCD42-10, scientific grade 1, MPP device		
	2048 x 512, 13.5 x 13.5 $\mu\text{m}$ pixels, 27.6 x 6.9 mm imaging area		
	<b>Typical</b>	<b>Maximum</b>	
<b>Dark Current @ -120°C</b>	0.3 e/p/hr	1 e/p/hr	
<b>System Read Noise</b>			
@ 100 kHz readout	3.5 e-rms	5 e-rms	
@ 1 MHz readout	8 e-rms	10 e-rms	
@ 2 MHz readout	13 e-rms	18 e-rms	
<b>Vertical shift rate</b> (software adjustable)	15.2 $\mu\text{sec}/\text{row}$		
<b>Spectral rate*</b>			
@ 100 kHz		35 spectra/sec	
@ 1 MHz		60 spectra/sec	
@ 2 MHz		90 spectra/sec	
	<b>Minimum</b>	<b>Typical</b>	
<b>Spectrometric Well Capacity</b>			
High Sensitivity	150 ke-	250 ke-	
High Capacity	600 ke-	800 ke-	
<b>Deepest Cooling Temperature</b>	-120°C	-110°C	
<b>Thermostat Precision</b>	$\pm 0.05^\circ\text{C}$ across entire temperature range		
<b>Software-selectable gains</b>	High	Mid	Low
High Sensitivity	1.5 e-/ct	3 e-/ct	6 e-/ct
High Capacity	6 e-/ct	12 e-/ct	24 e-/ct
<b>Dynamic Range</b>	16 bits		
<b>Nonlinearity</b>			
@ 100 kHz readout		< 1%	
@ 1 MHz readout		< 2%	
@ 2 MHz readout		< 2%	

\*Spectral rates measured with all rows vertically binned.

## QE Curves





www.piacton.com

email: moreinfo@piacton.com

USA +1.877.4 PIACTON | Benelux +31 (347) 324989

France +33 (1) 60.86.03.65 | Germany +49 (0) 89.660.779.3

UK +44 (0) 28.38310171 | Asia/Pacific +65.6293.3130

China +86 135 0122 8135 | Japan +81.3.5639.2741