

IsoPlane® 81

Aberration-Free Spectrograph
& Accessories

The Best Imaging Spectrograph on the Market



TELEDYNE PRINCETON INSTRUMENTS
Everywhere[™]youlook

Part of the Teledyne Imaging Group



Years of Research Satisfaction

When the first aberration-free IsoPlane 320 was introduced in 2013, it immediately became the **gold standard** for imaging spectrographs. Today, IsoPlane instruments remain the only aberration-free spectrographs on the market.

2013



"One of the ways we use IsoPlane 81 is together with a custom-made confocal microscope to collect single-particle fluorescence and scattering data that is then correlated with electron microscopy data."

—Dr. Denis Boudreau
Université Laval



"I can honestly say IsoPlane 81 is changing the way we do spectroscopy in our lab."

— Mark Waterland
Associate Professor
Massey University



"IsoPlane 81 offers a lot of convenience, all in one system"

—Dr. Nilam C. Shah
Lake Forest College



"The IsoPlane 81 system's small footprint and easy integration make it a perfect fit for our lab."

— Rajesh N. Davé,
Distinguished Professor
New Jersey Institute
of Technology

2020

Key Markets & Applications

Markets


Pharmaceutical

Drug discovery; PAT; QC;
Anti-counterfeit



Bio & Life Science

Cancer diagnostics; Surgery guidance;
Flow cytometry



Materials Science

2D; Nano; LED; Solar



Environmental Science

Droplet and aerosol; Pollutants;
Microplastics



Art & Archeology

Pigments; Conservation; Authenticity



Primary Applications

Raman
Spectroscopy
and Imaging


Fluorescence,
Emission,
Absorption

Hyperspectral
imaging

Microspectroscopy



#1 in Innovation for Highest SNR



ABERRATION-FREE SPECTROGRAPH

- No coma and astigmatism at all wavelength
- Diffraction-limited spectral resolution




DEEP COOLED CCD

- >95% QE
- Low dark current




eXcelon™

- Eliminates etaloning
- Improves QE



IntelliCal®

- Automated calibration
- NIST traceable ref. light sources
- Intensity and wavelength calibration



FAST DATA COLLECTION

- Dual readout ports
- Frame transfer
- Kinetics mode
- Mechanical shutter disabled mode

Perfect for OEM & System Integrators

Superior SNR with a footprint of a notebook. IsoPlane 81 is ideal to be used as an OEM Spectrometer where it can be easily integrated with other optical systems for Raman microscopy; hyperspectral imaging; cryo or elevated temperature measurement; flow cytometry; and more.

IsoPlane 81 comes with Labview®, Matlab®, and Python® drivers for easy software development and system integration.





IsoPlane 81 Models

	IsoPlane 81 BRX-VR	IsoPlane 81 BX-VR	IsoPlane 81 BRX-UR	IsoPlane 81 BX-UR
Sensor Type	Back-illuminated, deep-depletion, frame-transfer CCD with eXcelon™ technology and UV coating	Back-illuminated, frame-transfer CCD with eXcelon™ technology and UV coating	Back-illuminated, deep-depletion, frame-transfer CCD with eXcelon™ technology and UV coating	Back-illuminated, frame-transfer CCD with eXcelon™ technology and UV coating
Sensor Benefit	Enhanced NIR with >97% peak QE	Ultra-low dark current for long exposure experiment	Enhanced NIR with >97% peak QE	Ultra-low dark current for long exposure experiment
Spectral Range (nm)	400 - 1100	400 - 1100	200 - 1100	200 - 1100
Spectral Rate (burst mode)	>10,000 spectra/sec	>5,000 spectra/sec	>10,000 spectra/sec	>5,000 spectra/sec





Family Specifications

Focal Length	81 mm	
Aperture Ratio	f/4	
Spectral Resolution (FWHM)*	0.1 nm	
Spectrograph Optics & Spectral Range	VR: Protected silver coating 400 – 1100 nm	UR: UV-enhanced aluminum coating 200 – 1100 nm
Spatial Resolution	38.5 lp/mm @ 50% contrast over entire focal plane (Nyquist limited)	
Grating	150 g/mm up to 4320 g/mm. User changable, rotatable single-grating turret	
Astigmatism/Coma Aberration	Zero at all wavelengths and grating angles over entire focal plane	
Slits	Slit width: 10 μ m up to 500 μ m; 3.3 mm tall; interchangeable, laser-cut slits	
Camera	BRX	BX
Sensor Type	Back-illuminated, deep depletion, frame-transfer CCD sensor with eXcelon® and UV coating	Back-illuminated, frame-transfer CCD sensor with eXcelon® and UV coating
Sensor Format	1024 x 256 (1024 x 512 including frame-transfer storage area)	
Deepest Cooling Temperature	-55°C guaranteed, -60°C typical	
System Read Noise	4 e- @ 500 kHz; 6 e- rms @ 1 MHz; 27 e- @ 4.55 MHz	
Data Interface	USB 3.0 (3 m interface cable provided)	
Dimensions L x W x H	26.8 cm x 18.0 cm x 21.0 cm (11" x 7" x 8")	
Weight	8.84 kg (19.5 lbs)	

* at 750 nm with 10 μ m slit and a 1800 g/mm grating.

Specifications are subject to change.



>95%
QE



Ecosystem

IsoPlane 81 is much more than a spectrograph. Its meticulously engineered ecosystem is continually being expanded to support work in diverse fields. Fiber probes, accessories for Raman spectroscopy, and easy-to-use fiberoptics are just the beginning.



RAMAN MICROSCOPY

- CUBE-based standalone Raman microscope
- Compact footprint
- Built-in LED and witness camera



FIBER COUPLED PROBE

- Point-and-shoot
- Round-to-linear collection fiber bundle for best efficiency
- Fiber adaptor can accept different fiber input including customized probe



RAMAN MICROSCOPY

- Convert standard microscope to Raman microscope
- Plug-and-play
- Capable of confocal measurement



CUVETTE HOLDER CONFIGURATION

- Ideal for routine sample measurement
- Accepts standard 10 mm cuvette
- Maintenance-free



There's a CUBE for That!

CUBES



**Lasers and
calibration
light sources**



Designed
for free
space and
fiber
coupling

Pre-aligned
(stays aligned
as more
CUBES are
added)

With IsoPlane 81 CUBES & accessories, interfacing is a cinch!

If you've ever spent more time digging through cabinets than conducting your actual experiment, you'll definitely appreciate these CUBES.



Easy to
attach using
only a
screwdriver

Compatible
with
Thorlabs® 30
mm cage
system

**User
changeable
slits and
grating**








And more...





LightField Enables Better Outcomes

While IsoPlane 81 successfully distinguishes itself with innovative hardware, the beautifully designed LightField software orchestrates all inner workings and makes it a pleasure to use the system day in and day out!





Dedication

Customer for Life

IsoPlane 81 users look at biological or physical processes using optical spectroscopy techniques. They may work at a university, a private company, or a government laboratory. They may be a scientist, an educator, or an engineer. Their level of experience as a spectroscopist may vary, too.

IsoPlane 81 ecosystem allows them to perform advanced spectroscopic research without needing to integrate countless third-party components. The aberration-free results yielded by the system's fully integrated hardware and software are worthy of peer review.



Quality

Teledyne Princeton Instruments builds quality and reliability into every one of our products, including a lifetime warranty on vacuum systems.

Partnership

From Nobel laureates to Original Equipment Manufacturers (OEM), we dedicate ourselves to understanding the unique requirements of every customer.

Knowledge

Teledyne Princeton Instruments employs a diverse and knowledgeable staff from all disciplines and backgrounds.

Support

Our worldwide locations and experienced technicians support you before, during, and especially after the sale.





Who We Are

3000+


employees worldwide




TELEDYNE TECHNOLOGIES
Everywhere you look™

When you partner with Teledyne
Princeton Instruments, you become part
of a world-class organization with
diversified technologies and capabilities.

Machine
Vision



Scientific
Imaging




MEMS




Healthcare



Astronomy &
Earth



Geospatial



www.princetoninstruments.com
pi.info@Teledyne.com
+1 609.587.9797

For information regarding any of the images/data shown in this brochure, please contact Teledyne Princeton Instruments.

Copyright © 2020 Teledyne Princeton Instruments. All rights reserved. IsoPlane, LightField, and NIRvana are registered trademarks of Teledyne Princeton Instruments. Linux is the registered trademark of Linus Torvalds in the U.S. and other countries. Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and other countries. Python is a registered trademark of the Python Software Foundation. All other brand and product names are the trademarks or registered trademarks of their respective owners and manufacturers.

Rev0



TELEDYNE PRINCETON INSTRUMENTS
Everywhereyoulook™

Part of the Teledyne Imaging Group