

High-performance near-infrared Raman for clinical application

PENG ZOU

Technology advances are enabling Raman spectroscopy for application as a clinical tool.

Raman spectroscopy has seen tremendous growth in biological and life science applications over the past two decades. As it measures inelastic scattering from a monochromatic light source incident onto a sample, Raman spectroscopy provides rich information on molecular structures, identities, and composition, and can be used for both qualitative and quantitative chemical analysis. This information, indicated by the wavelength difference (or relative Raman shift) and intensity of the scattering, results from interactions between photons and molecular vibrations.

This label-free, noninvasive technique has been widely used in research, and has more recently been pursued for its potential as an alternative or complementary clinical tool for patient care and disease diagnostics. This article will discuss the unique requirements for a medical Raman instrument, recent innovations in hardware development, and clinical cancer diagnostic applications enabled by these technologies.

Raman spectroscopic systems

Compared to other optical spectroscopy methods, Raman spectroscopy is considered a weak phenomenon with a

relatively low cross-section ($<1 \times 10^{-8}$ ratio of Rayleigh vs. Raman scattering). High-sensitivity Raman spectrometers have been at the heart of successful research and application development. A Raman spectrometer comprises a laser source, spectrograph, and a detector (see Fig. 1). Although free-space coupling is used in many Raman systems, fiber-optic probes are widely used in medical diagnostics because of their flexibility and collection efficiency enabled by fiber bundling. Progress in hardware and software development—including high-throughput spectrographs, sensitive

with high accuracy and high sample throughput at a much lower cost with no sample preparation—capabilities that are either technically impossible or time- and cost-prohibitive using conventional methods (see Table 1). With the ability to monitor various diseases and key health indicators, Raman spectroscopy can provide real-time information for surgical guidance, drug effectiveness measurement, and point-of-care diagnostics when aided by an appropriate novel sampling interface.

Designing a Raman spectrometer for a clinical environment requires consideration of the unique optical properties of biological tissues:

- Very high water content
- Highly fluorescent in the visible (VIS) region
- Highly scattering
- Strongly absorbent in the UV-VIS region

Excitation laser

Autofluorescence from biological specimens introduces severe interference to Raman spectral data and obscures analysis of results. However, longer-wavelength excitation lasers can effectively reduce or eliminate it. As longer wavelengths penetrate much deeper than UV-VIS light, Raman spectral data collected with near-infrared (near-IR) lasers have more biochemical information needed for accurate *in vivo* under-skin, tissue, and tumor analysis. Two optical windows in the near-IR spectral range are recommended




FIGURE 1. The key components of a Raman spectroscopy system are the laser, detector, spectrograph, and probe.

detectors with high quantum efficiency and low dark noise, novel sampling interfaces, and robust chemometrics tools—have made Raman spectroscopy an increasingly practical and powerful tool for complex biological cell- and tissue-related research.

Compared to conventional diagnostic methods used in medicine, Raman spectroscopy has potential to enable both *in vivo* and *ex vivo* diagnostics

Table 1. Comparison of Raman spectroscopy with conventional disease diagnostic methods

	Raman spectroscopy	Conventional
Technology	Laser-based	Biopsy; endoscopy; MRI; CT; x-ray
Throughput	Fast	Time-consuming
Cost	Low	High
Speed	Real time	Slow
Sampling	Noninvasive	Invasive
Regulatory	Developing	Approved

for Raman spectroscopy (see Fig. 2). The first, and most commonly used, lasers for Raman excitation are 785 or 830 nm lasers, which overlay well within the first optical window. The second are 1064 nm lasers, which are used when autofluorescence is severe with shorter wavelengths. The choice of laser also relies heavily on the appropriate detector, which will be discussed in the next section.

For *in vivo* measurements, the selection of laser power is determined by a few factors:

- Maximum permissible exposure
- Best signal intensity within the safety allowable power level
- Comfort level
- Local heating effect

For *ex vivo* measurements, a higher laser power may be used to improve the signal intensity.

Since biological tissues can have complex textures and can be highly scattering, both multi- and single-mode lasers may be used for a clinical Raman system. The latter are preferred for Raman microscopic

Optical penetration depth (mm)




FIGURE 2. Optical windows of Raman spectroscopy in biological tissue. (Adapted from A. N. Bashkatov [5])




FIGURE 3. The near-IR optimized TPIR 785 Raman spectrometer (Teledyne Princeton Instruments) achieves better light throughput than a commercial lens spectrometer, as demonstrated in this comparison using the 1033 cm⁻¹ Raman peak of methanol with a 785 nm laser.

measurements when high spatial resolution is required.

While near-IR light greatly reduces fluorescence interference, signal strength decreases with longer wavelength excitation because the spontaneous Raman scattering cross-section is proportional to $1/\lambda^4$, which makes the already ‘weak’ Raman signal even weaker in the near-IR. Also, water has a small Raman cross-section and human tissue is rich with water. This further reduces the Raman scattering from human tissue. One can potentially increase laser power to improve Raman signal to certain level, but safety concerns and highly photon-sensitive tissues prohibit using a high-power laser. All these challenges point to a high-performance

near-IR Raman spectrometer as the fundamental requirement for a successful clinical instrument.

Spectrograph

Dispersive grating-based spectrographs have several advantages over Fourier transform interferogram-based systems, and they are more broadly used as medical diagnostic tools due to simple optical design, lack of moving parts, and fast data collection with a focal-plane array detector.

Biological tissues are composed of mainly water with large molecules of proteins, lipids, and carbohydrates. Therefore, the Raman spectra from biological samples typically have broad peaks and features and do not need a long focal-length and high-resolution spectrograph.

On the other hand, a spectrograph with high throughput of light and low f/N is more desirable for high signal-collection efficiency. In addition, high-quality antireflective coating of the optics, aberration-corrected spectrograph design, and high efficiency gratings are critical to achieve the highest light throughput.

As an example, Teledyne Princeton Instruments’ TPIR-785 Raman spectrometer incorporates a near-IR optimized lens spectrograph and thus easily enables a 40–60% light throughput improvement throughout the spectral range (see Fig. 3).

The TPIR-785’s spectrograph uses custom designed lens optics for aberration correction and provides imaging quality far superior to that of conventional Czerny-Turner spectrographs (see Fig. 4).




FIGURE 4. Lens optics designed for aberration correction enable improved imaging quality; images of a 19 × 200 μm fiber array of 880 nm line (top) and an image of a 19 × 200 μm fiber array with a broadband light source (bottom) are shown.

THE STEEPEST
EDGE FILTERS
ON THE MARKET
DOWN TO 20 cm^{-1}


The Leader in Optical Filters for Raman Applications

THREE FAMILIES OF PRODUCTS TO BEST MEET YOUR PERFORMANCE AND COST NEEDS

- › View the weakest signals closer to the laser line than you ever have before
- › With their deep laser-line blocking, ultra-wide and low-ripple passbands, proven hard-coating reliability, and high laser damage threshold, our filters offer performance that lasts
- › Edge wavelengths available from 224 to 1550 nm
- › Adjust the angle of a Raman filter to shift its edge wavelength for customized performance – model the shift at <https://searchlight.semrock.com>

Learn more at www.semrock.com/raman-spectroscopy.aspx

Edge Filter Type	Guaranteed Transition Width (% of laser wavelength)	Typical Edge Steepness (% of laser wavelength)
RazorEdge E-grade	< 0.5% (63 cm^{-1} for 785)	< 0.2% (25.4 cm^{-1} for 785)
RazorEdge U-grade	< 1.0% (126 cm^{-1} for 785)	< 0.5% (63 cm^{-1} for 785)
EdgeBasic	< 2.5% (309 cm^{-1} for 785)	< 1.5% (186 cm^{-1} for 785)



Transition Width = maximum allowed spectral width between the laser line (where OD > 6) and the 50% transmission point. Edge Steepness = actual steepness of a filter measured from the OD 6 point to the 50% transmission point.

Direct benefits of such high-quality imaging include:

- Multichannel or Raman imaging measurement with minimum crosstalk
- Stable high etendue and resolution across the entire spectral range and detector focal plane
- Works seamlessly with a Raman probe with a large fiber bundle and truly benefits from state-of-the-art large-format detectors

Detectors and probes

For high-performance spectrometers, charge-coupled device (CCD) detectors still provide the best quantum efficiency (QE) and lowest dark current with proper cooling in the UV-VIS region (see Fig. 5). However, the CCD QE starts to decrease at ~800 nm and has a cutoff at ~1100 nm, which is not quite efficient for near-IR Raman measurements. A high QE detector is desired for 785–1100 nm in order to effectively cover the fingerprint Raman spectral region (150–1500 cm^{-1}). The TPIR-785 has the option to select a recently developed super-deep-depletion CCD camera (BLAZE HR) as its detector. The depletion region of epitaxial silicon is much thicker than the typical back-thinned CCD that drives QE in the near-IR region to 2–7X higher than conventional CCD detectors. The state-of-the-art camera




FIGURE 5. Quantum efficiency (QE) curves of a conventional back-thinned deep-depletion CCD (blue); super-deep-depletion CCD (BLAZE HR camera from Teledyne Princeton Instruments, orange); and InGaAs detectors (black) are shown.

also provides deep cooling to -95°C with air for lowest dark current. The high QE in near-IR region and low dark current makes BLAZE HR optimal for near-IR Raman. Additionally, the camera provides ultra-high readout speed through dual ports for fast spectral rate, which is critical for real-time diagnostics and Raman imaging.

Some research areas, using Raman spectroscopy, have utilized 1064 nm lasers because they overlay with the second optical window and provide even less fluorescence background. Indium gallium arsenide (InGaAs) sensors are typically paired with 1064 nm

lasers to provide better coverage in longer-wavelength (900–2200 nm) regions. But InGaAs detectors suffer from greater dark noise than silicon-based CCD detectors, and with a further Raman cross-section decrease at 1064 nm, as well as eye-safety concerns due the invisibility of the laser






FIGURE 6. Raman probes for clinical diagnostics include one with a small probe tip



that enables precise location information for in vivo diagnostics and boundary detection during a surgery (a) and a multimodality probe head that allows simultaneous Raman and diffuse reflective measurements for more accurate diagnostic results (b). (Courtesy of EmVision LLC)



CRISP Autofocus System



The CRISP substantially eliminates focus drift in high-power microscopy applications by sensing minute changes between the objective lens and the specimen cover slip. It also allows a specimen to remain accurately focused for hours at a time with an accuracy of 5% of the objective depth of focus and maintains focus while scanning in XY.

www.asiimaging.com
info@asiimaging.com
541.461.8181
800.706-2284

beam, 1064 nm Raman systems are not widely used for medical diagnostics.

Another critical component for clinical Raman systems is the optical fiber probe. Raman probes not only provide flexible operation and easy access to different parts of the body for *in vivo* measurements, they also can be used to improve Raman signal collection efficiency and expand measurement capability. Such novel probes enable endoscopic Raman, Raman imaging, and multimodality measurement (see Fig. 6).^{1,2}

Cancer diagnostics and other applications

The noninvasive, simple, fast, and highly reproducible characteristics of its measurements make Raman spectroscopy a strong candidate for high-profile applications including cancer diagnostics, microbiology identification, and noninvasive blood glucose measurement.^{3,4}

Table 2. Overview of Raman-based techniques

	Advantages	Disadvantage
Spontaneous Raman spectroscopy	Most broadly studied; simple hardware setup; no sample preparation	Weak signal; requires high-performance instrument
Surface-enhanced Raman spectroscopy (SERS)	>10 ¹⁰ Raman cross-section enhancement; greatly improves the sensitivity	Development of biocompatible SERS substrates, which are typically made of metal nanoparticles; stability and repeatability
Coherent anti-Stokes Raman spectroscopy (CARS)	Orders of magnitude larger cross-section; Excellent for imaging	Complex and costly system setup involving multiple lasers
Spatially offset Raman spectroscopy (SORS)	Provides depth resolution into tissues by collecting scattering at different distances from the laser incident point	Requires high-sensitivity Raman system; the unique SORS fiber bundle probe also requires large-format camera with tall slit
Resonance Raman spectroscopy	Greatly enhanced intensity of the Raman scattering	Needs to meet the enhancement condition; not for all samples or chemicals

NEW! Open Design Upright Microscope

- Simple coarse height adjustment
- Built-in motorized focus
- *In vivo* and *In vitro* experimentation on one setup
- OCC or DIC transmitted light (LED)
- Epi-fluorescent imaging
- Uses standard Olympus objectives



BOB

The Sutter BOB – designed to eliminate the conventional microscope frame – is a versatile, open-design upright microscope platform ideal for slice electrophysiology, widefield fluorescent imaging, two-photon imaging, photostimulation and new techniques just being developed!

SUTTER INSTRUMENT

PHONE: +1.415.883.0128 | FAX: +1.415.883.0572
EMAIL: INFO@SUTTER.COM | WWW.SUTTER.COM

A large amount of research, covering nearly every aspect of medical diagnostics, has been dedicated to clinical Raman spectroscopic applications over the past decade. Pioneering studies using Raman spectroscopy for detection of cancerous tissue can be found as early as the 1990s.^{6,7} Since then, cancer diagnostics has become one of the most prolific research areas:^{4,2}

The number of publications on Raman spectroscopy for cancer diagnostics has grown from less than 4000 in 2010 to more than 14,000 in 2019. Highlights in this field include:

- Diagnostics using Raman spectroscopic methods performed on different types of cancer; some results were very promising for early detection

- *In vivo* measurements showing potential for real-time diagnostics
- *Ex vivo* disease detection developed based on Raman spectroscopy from biopsy samples, body fluid, and blood
- Clinical trials
- The development of novel chemometrics, machine learning, and artificial intelligence methods to build robust prediction models

Technology development greatly improves the feasibility and viability of Raman spectroscopy for clinical applications. Since the 1928 discovery of the Raman phenomenon, many innovations have been developed to improve certain aspects of the technique (see Table 2). Improvement of instrument performance will benefit all Raman-based methods.

More is needed, however, for Raman to fulfill its promise. It is interesting to reflect that the modernization of pharmaceutical manufacturing, through the adoption of optical spectroscopic process analysis tools required a joint effort involving research communities, equipment vendors, pharmaceutical companies, and regulatory bodies for decades before the technologies finally became viable for the industry.

Among near-term challenges that will help Raman spectroscopy gain acceptance for clinical application are cost reduction; robust chemometrics/artificial intelligence/machine learning tools; size, weight, and power reduction; and regulatory acceptance. ◀

REFERENCES

1. K. Lin et al., *Theranostics*, 7, 14, 3517–3526 (2017).
2. E. Cordero et al., *J. Biomed. Opt.*, 23, 7, 1-23 (2018).
3. H. J. Butler et al., *Nat. Protoc.*, 11, 664–687 (2016).
4. I. Pence et al., *Chem. Soc. Rev.*, 45, 1958–1979 (2016).
5. A. N. Bashkatov, *J. Phys. D: Appl. Phys.*, 38, 2543 (2005).
6. M. S. Feld et al., *Advances in Fluorescence Sensing Technology II*, 2388, 99–104 (1995).
7. A. Mahadevan-Jansen et al., *International Conference of the IEEE Engineering in Medicine and Biology Society*, 6, 2722–2728 (1997).

Peng Zou, Ph.D., is Product Manager at Teledyne Princeton Instruments, Trenton, NJ; e-mail: peng.zou@teledyne.com; princetoninstruments.com.



BWTEK
A Metrohm Group Company

Portable Raman Analyzer for Measurements Through Opaque Packaging

The **STRam**® provides easy identification of materials through a variety of packaging and barrier layers utilizing 785nm or 1064nm wavelengths!

Learn More about the STRam
www.bwtek.com/STRam

+1-302-368-7824 www.bwtek.com marketing@bwtek.com