



## NIRvana® 640

#### **SWIR For Science**

- Unbeaten Signal to Noise Ratio
- Ultra-Deep Cooling for Long Exposures
- Lowest Dark Current Noise Without Liquid Nitrogen
- Cold Shielding to Minimise Scene Noise



#### What Makes Scientific Quality SWIR Cameras?

With Scientific Imaging comes unique challenges. Light levels are often low, image quality requirements are high, and reliability is paramount.

The NIRvana® range of Scientific InGaAs cameras from Teledyne Princeton Instruments is specifically designed to meet the demanding requirements of scientific imaging.


With sensitivity often the primary requirement for scientific imaging, the NIRvana range uses large pixels, high quantum efficiencies and unbeaten low dark current noise to provide industry-leading signal to noise ratios. Combined with this is the incredible image quality achieved through advanced pixel correction algorithms, overcoming the inherent challenges of InGaAs sensors.

Camera noise sources are not alone in infrared imaging. Thermal photons from the scene and from the camera housing are everpresent, conveying background noise. The NIRvana series uses integrated cold shielding to keep these photons, and hence this unwanted background noise out of images.


High speeds and short minimum exposure times provide the versatility required for scientific imaging. From whole animal imaging to spectroscopy of nanomaterials, the NIRvana range has a solution for your application.

#### Further Benefits of the NIRvana 640:


- High Dynamic Range with >600,000e- Full Well Capacity
- Air cooling for simplicity, or liquid cooling for vibrationsensitive applications
- Lifetime Vacuum Guarantee
- Imaging and spectroscopy possible within one powerful, intuitive software package




Sunspot Imaging



Imaging &
Spectroscopy of
Nanomaterials



In vivo Imaging



Singlet Oxygen Imaging



#### NIRvana Application Highlights

NIRvana® has a peerless scientific pedigree with **10 publications in Nature journals using NIRvana® cameras since 2017**. Whatever your scientific interest, and whether your research takes place in a university lab or in industry, the NIRvana range of cameras is there to deliver the performance you need.

NIRvana cameras cover a massive range of application areas. Below are some comments shared with us by researchers using NIRvana in their experiments.


- Nanotube fluorescence
- Small-animal imaging
- Nondestructive testing
- Semiconductor analysis
- NIR probe development
- NIR-II imaging
- NIR spectroscopy
- 1064 nm Raman spectroscopy
- Astronomy (J band, H band)
- Adaptive optics
- Low-light / night vision
- Laser-illumination imaging

#### Deep Penetration in vivo Imaging in NIR-II

"The NIRvana 640 is highly sensitive, with a high signal-to-noise ratio, providing clear structural information about the blood vessels."

Prof. Fan Zhang, & Dr. Zhang Hongxin, Department of Chemistry, Fudan University, China





#### Measuring Spectra of Single Quantum Dots

"The unprecedented long integration times and low dark counts [of the NIRvana LN] enabled single-nanocrystal spectral experiments, which cannot be achieved with conventional InGaAs detectors for these dots."


Dr. Han Htoon, Center for Integrated Nanotechnologies, Los Alamos National Laboratory

#### **Energy-Momentum Resolved Spectroscopy**

"Unlike other InGaAs 2D detectors we tested, the TE cooling and background subtraction features have allowed us to easily acquire long exposures to resolve photoluminescence from very dilute emitter samples."

Prof. Rashid Zia,

School of Engineering, Brown University, USA






### What Sets NIRvana Apart

## Lowest Dark Current Without Liquid Nitrogen


The NIRvana 640's multi-stage thermo-electric cooling (TEC) reaches -85°C with water cooling, reducing dark current noise to only 40e-/p/s – two order of magnitude lower than industrial cooled SWIR cameras.

As a result, the sensitivity of the NIRvana 640 is absolutely without match for air or water cooled SWIR cameras.

The Liquid Nitrogen-cooled NIRvana LN is the only camera that can rival the NIRvana 640's sensitivity.



Single-Walled Carbon Nanotube (SWCNT)<sup>2</sup>



## Unparalleled Thermal Management

The exemplary thermal management and lifetime vacuum guarantee of NIRvana cameras demonstrate Teledyne Princeton Instruments' 40 years of experience in cooling technology.

The low temperatures achieved, and as importantly the stability with which they are maintained, provide the lowest noise for long exposures and vital baseline stability for quantitative measurements<sup>1</sup>.

Our all-metal hermetic vacuum seals are completely maintenance-free. The NIRvana series are the only SWIR cameras guaranteed against moisture or loss of vacuum for the entire lifetime of the camera.

<sup>&</sup>lt;sup>2</sup> - Image courtesy of Dr. R. Bruce Weisman, Professor of Chemistry, The Richard E. Smalley Institute for Nanoscale Science and Technology, Rice University, Houston, TX USA





<sup>&</sup>lt;sup>1</sup> - For further information, please see our NIRvana Thermal Management Technical Note: www.princetoninstruments.com/nirvana/thermal


#### **Optimum Image Quality**

InGaAs sensors are subject to a number of factors that affect image quality. Individual pixel defects, star defects, column-to-column variations and multiple pixel clusters are all universal features of InGaAs sensors. Without careful correction, these defects can throw off analysis routines and reduce image quality.

The NIRvana range of cameras have unique capabilities for dealing with these defects powerfully and transparently<sup>3</sup>, providing a scientific image quality and imaging workflow that provides data you can trust.







# Unwanted SWIR Photons from Scene Photons from Scene

Sensor

Cold Shield

**SWIR Photons** 

from Camera Electronics

#### Minimized Scene Noise

A unique challenge in infrared imaging is that SWIR photons are, in fact, everywhere. We aren't of course able to see them to allow us to easily trace and eliminate background photon sources. But all background light contributes background noise, which can in many cases outstrip camera-based noise sources depending on the severity. Indeed, due to blackbody radiation, all room temperature objects are constantly emitting SWIR photons, as are the electronic circuits and housings of SWIR cameras.

The NIRvana range of cameras features Cold Shielding. This is a physical barrier that blocks unwanted photons outside a specified cone angle,  $\theta$ , as shown in the diagram opposite. This shield is cooled to the same ultra-low temperature as the sensor leading to negligible thermal photon emission<sup>4</sup>.

Additionally, world-leading thermal management engineering isolates thermal photons from camera electronics and maintains a lower camera housing temperature, reducing further the possibility of unwanted photon detection.

<sup>&</sup>lt;sup>4</sup> - For further information, please see our NIRvana Cold Shielding Technical Note: <u>www.princetoninstruments.com/nirvana/shielding</u>




<sup>&</sup>lt;sup>3</sup> - For further information, please see our NIRvana Image Corrections Technical Note: <u>www.princetoninstruments.com/nirvana/corrections</u>

## **Unbeaten Sensitivity**

The NIRvana range of cameras represent the height of SWIR sensitivity.

#### Large Pixels Collect More Light

As sensitivity scales with pixel area, the 20  $\mu$ m pixels used by the NIRvana range of cameras provide 1.8x the light collection ability of a 15  $\mu$ m pixel, and 4x that of a 10  $\mu$ m pixel.




#### Low Noise Operation

Dark Current Noise is typically the largest noise source in InGaAs imaging. Due to the low cooling temperatures achieved, dark current noise is drastically reduced in the NIRvana range of cameras. Combined with low read noise, the NIRvana range provides incredible signal to noise performance.

#### **High Quantum Efficiency**

The NIRvana 640 achieves >80% quantum efficiency across an incredibly broad range, from 950 to 1600nm, enabling low light imaging across the NIR-II / SWIR region.





#### Combining These Factors: Signal to Noise Ratio

What is the relative importance of pixel size, quantum efficiency, dark current, read noise, and all the other factors above? The answer depends upon your typical exposure time, background light level, and expected photon budget. But all the complexity can be set aside: **the Signal to Noise Ratio represents the entire balance of factors with one single number**<sup>5</sup>. Equation 1 shows how this is calculated:

$$SNR = \frac{Signal}{Total\ Noise} = \frac{I_s \times t_{exp} \times d_{pixel}^2 \times QE}{\sqrt{\left((I_s + I_b) \times t_{exp} \times d_{pixel}^2 \times QE\right) + \left(\sigma_d \times t_{exp}\right) + \sigma_r^2}}$$

Equation 1: Signal to Noise Ratio for InGaAs Cameras. The first term in Total Noise represents total photon shot noise.

 $I_s$  = Irradiance, i.e. incoming signal in (photons / second /  $\mu$ m<sup>2</sup>)

 $t_{exp}$  = Exposure time (seconds)

 $d_{pixel}$  = Pixel size (µm)

QE = Quantum Efficiency at the detection wavelength of interest

 $I_b$  = Background irradiance, i.e. all background signal (photons / second /  $\mu$ m<sup>2</sup>)

 $\sigma_d$  = Dark current (e- / pixel / second)

 $\sigma_r$  = Read noise (e-)

The graphs on the following page show the calculated signal to noise ratios of the three NIRvana cameras compared to three representative alternative cameras. For the calculation the contribution of scene noise from the experimental setup,  $I_b$ , is ignored.


Incoming signal is explored over 3 orders of magnitude from 'high' signals of 10 photons / s /  $\mu m^2$ , which would equate to 4000 photons per pixel per second for the NIRvana cameras, to 'very low' signals of 0.01 photons / s /  $\mu m^2$ , for which very long exposure times would be necessary to achieve adequate signal to noise. The graphs are shown up to the point where pixel saturation would occur through the combination of collected signal and dark current.

The dominance of low dark current and large pixel sizes is evident; **the NIRvana camera range offers superior sensitivity in practically every case.** For the alternative cameras, due to their higher dark current, saturation of the pixel is frequently seen before an adequate signal to noise ratio can be achieved.

In real-life imaging situations, scene noise that we discounted here will be present in almost all cases. The presence of **cold shielding** in the NIRvana range of cameras will then be very likely to **significantly increase the NIRvana cameras' lead in signal to noise ratio**, compared to alternatives without cold shielding which would suffer from greater background noise  $I_b$ .



<sup>&</sup>lt;sup>5</sup> - For further information, please see our NIRvana Signal and Noise Technical Note: www.princetoninstruments.com/nirvana/signalandnoise



Signal to Noise Ratio graphs for different incident signal levels, defined in photon flux, in photons per second per square micrometer incident on the pixel. Each SNR line ends at pixel saturation. The NIRvana range of cameras are shown in comparison to three alternative cameras representative of cameras on the market, with specifications given in the table above. Theoretical calculation based on specification sheet values.



## NIRvana 640 Specifications

Feature	NIRvana 640	
Sensor	2D InGaAs FPA	
Imaging format	640 x 512	
Active area	12.8 mm x 10.24 mm	
Pixel size	20 μm x 20 μm	
Optical Mount	C-Mount, Spectrometer Mount	
Spectral response	0.9 μm to 1.7 μm	
Quantum efficiency	>80% from 1µm to 1.6µm; peak QE is 85%	
Pixel well capacity (e-)	High Gain ≥ 45,000	High Capacity ≥ 600,000
Nominal gain (e-/count)	1	14
Typical system read noise (e-)	<75 (High Gain)	
Response nonlinearity <sup>6</sup>	<2% (High Gain)	
Cooling temp. @ 20°C ambient	-80°C air (guarranteed), -85°C (typica	nl) with Air + 15°C liquid circulation
Cooling method	Air only, liquid only, or a combination of air and liquid	
Dark current (e-/p/s) <sup>7</sup>	40 @ -85°C, 150 @ -80°C	
Thermostating precision	±0.05°C across entire temperature range	
Blemish specification <sup>8</sup>	Grade A: <2% defects	
Cold shield <sup>9</sup>	f#/1.5; Light cone angle $\theta = 20^{\circ}$	
Output interface	GigE	
Digitization	16 bits	
Scan rate	10MHz, 5MHz, 2MHz	
Frame rate	110 fps @ 10MHz 55 fps @ 5MHz 22 fps @ 2MHz	
Binning and ROI	Software only	
Shutter	Global shutter	
Exposure time	<2 µs to 10 minute	
Window material	Fused silica (AR coated)	
Operating temperature	0°C to +30°C	
Weight	9.5 lbs (4.3 kg)	
Dimensions (C-mount) (L x W x H)	7.40" x 5.8" x 5.8" (187.96 mm x 147.32 mm x 147.32 mm)	
Certification	CE	

<sup>&</sup>lt;sup>6</sup> - Linearity for exposure > 20 ms

All specifications are typical performance and are subject to change. NOTE: Export of this camera outside of the United States is prohibited by law unless accompanied by a valid Export License as issued by the United States Department of Commerce.




<sup>&</sup>lt;sup>7</sup> - Measured with a cold target at -174°C (99K)

<sup>&</sup>lt;sup>8</sup> - For detailed blemish specifications, please contact Princeton Instruments

 $<sup>^{\</sup>rm 9}$  - See page 5 for explanation of Cold Shield specifications

#### Software Support & Integration


LightField® Software allows for complete control of all Teledyne Princeton Instruments cameras and spectrometers on an easy to use platform. LightField software also includes simple image post-processing software and a built-in, smart math engine to obtain the most from acquired data.



Teledyne Princeton Instruments LightField® provides imaging and spectroscopy within one powerful and versatile software package.

LightField® also makes integration easy with LabVIEW® (National Instruments), MATLAB® (MathWorks) or Python®, allowing you to command the software directly from whichever program you prefer.

For those requiring a custom solution, the PICAM API and SDK, from Teledyne Princeton Instruments offer complete control over all cameras, available for all 64-bit Windows and Linux systems.





#### Teledyne Princeton Instruments: Part of Teledyne Imaging

Teledyne Imaging is a group of leading-edge companies aligned under the Teledyne umbrella. Teledyne Imaging forms an unrivaled collective of expertise across the spectrum and decades of experience. Individually, each company offers best-in-class solutions. Together, they combine and leverage each other's strengths to provide the deepest, widest imaging and related technology portfolio in the world.


Teledyne Princeton Instruments (TPI), a business unit of Teledyne Digital Imaging US, Inc, designs and manufactures high-performance CCD, CMOS, ICCD, EMCCD, emICCD, and InGaAs cameras; spectrographs; and optics-based solutions for the scientific research, industrial imaging, and OEM communities. We take pride in partnering with our customers to solve their most challenging problems in unique, innovative ways.

Teledyne Princeton Instruments has a vast history of facilitating key research in many fields. Our family of NIRvana® cameras has enabled numerous scientific and commercial advances by delivering the world's best SWIR imaging performance.


#### The NIRvana Camera Range

Need even more speed? Consider the NIRvana HS!

Need even more sensitivity? Consider the NIRvana LN!




#### Outline Drawings with C-Mount



Spectrometer adapter is available upon request.

#### **Outline Drawings with Spectrometer Mount**



C-mount adapter is available upon request.



## NIRvana® 640 SWIR FOR SCIENCE



Contact your local Teledyne Princeton Instruments representative for additional information.

#### Teledyne Princeton Instruments – USA

Tel: +1 609.587.9797 pi.info@teledyne.com

#### China

Tel: +86 157 2153 5343 pi.info.china@teledyne.com

#### France

Tel: +33.1.70.38.19.00 EVR@teledyne.com

#### Germany

Tel: + 49 (0) 89-660 779 3 pi.germany@teledyne.com

#### Japan

Tel: +81.3.6709.0631 Pl.Nippon@teledyne.com

#### **United Kingdom**

Tel: +44 (0) 7810 835 719 pi.info@teledyne.com

Copyright® 2021 Princeton Instruments, Inc. All rights reserved. NIRvana, eXcelon, and LightField are registered trademarks of Princeton Instruments, Inc. Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and other countries. Linux is the registered trademark of Linus Torvalds in the U.S. and other countries. Python is a registered trademark of the Python Software Foundation. All other brand and product names are the trademarks or registered trademarks of their respective owners and manufacturers.

Use and Disclosure of Data Information contained herein is classified as EAR99 under the U.S. Export Administration Regulations. Export, re-export or diversion contrary to U.S. law is prohibited.



