

# Lansis® System Manual



4411-0174  
Issue 2  
July 19, 2022

[www.princetoninstruments.com](http://www.princetoninstruments.com)

Use and Disclosure of Data Information contained herein is classified as EAR99 under the U.S. Export Administration Regulations. Export, reexport or diversion contrary to U.S. law is prohibited.

## Revision History

Issue	Date	List of Changes
1	July 19, 2022	<ul style="list-style-type: none"><li>Initial version of this document.</li></ul>

©Copyright 2017-2022    Teledyne Princeton Instruments  
All Rights Reserved    3660 Quakerbridge Rd  
                                 Trenton, NJ 08619  
                                 TEL: 800-874-9789 / 609-587-9797

All rights reserved. No part of this publication may be reproduced by any means without the written permission of Teledyne Princeton Instruments.

Printed in the United States of America.

Lansis, IntelliCal, IsoPlane, LightField, and PVCAM are registered trademarks of Teledyne Digital Imaging US, Inc..

LabVIEW is a registered trademark of National Instruments, Inc.

LEMO is a registered trademark of INTERLEMO HOLDING SA

Windows is a registered trademark of Microsoft Corporation in the United States and/or other countries.

The information in this publication is believed to be accurate as of the publication release date. However, Teledyne Princeton Instruments does not assume any responsibility for any consequences including any damages resulting from the use thereof. The information contained herein is subject to change without notice. Revision of this publication may be issued to incorporate such change.

# Table of Contents

---

<b>Chapter 1:</b>	<b>About this Document.....</b>	<b>7</b>
1.1	Intended Audience.....	7
1.2	Related Documentation.....	7
1.3	Document Organization.....	8
1.4	Safety Related Symbols Used in this Manual.....	9
1.5	Lansis Safety Information.....	10
1.6	Precautions.....	11
<b>Chapter 2:</b>	<b>Lansis Camera System.....</b>	<b>13</b>
2.1	Lansis Camera.....	14
2.1.1	eXcelon®.....	14
2.1.2	Power.....	15
2.1.3	CCD Arrays.....	16
2.1.4	Cooling.....	16
2.1.4.1	Internal Fan.....	16
2.1.4.2	External Cooling Circulator.....	16
2.1.4.3	Coolant Ports.....	17
2.1.5	Rear-Panel Connectors and Indicators.....	18
2.2	Cables.....	19
2.3	Certificate of Performance.....	19
2.4	Application Software.....	20
2.5	Accessories.....	21
2.5.1	External Shutters.....	21
2.6	Unpack the System.....	21
2.6.1	Verify Equipment and Parts Inventory.....	21
2.7	Lansis Camera and System Maintenance.....	22
2.7.1	Camera.....	22
2.7.2	Optical Surfaces.....	22
2.7.3	Repairs.....	22
<b>Chapter 3:</b>	<b>Install LightField.....</b>	<b>23</b>
3.1	Prerequisites.....	23
3.2	Installation Procedure.....	23
<b>Chapter 4:</b>	<b>System Block Diagrams.....</b>	<b>25</b>
<b>Chapter 5:</b>	<b>Hardware Configuration.....</b>	<b>27</b>
5.1	Spectrograph Adapters.....	27
5.2	Overexposure Protection.....	27
5.3	External Chiller Use.....	28
5.3.1	Experiment Shutdown.....	29
<b>Chapter 6:</b>	<b>LightField First Light.....</b>	<b>31</b>
6.1	Set Up and Configuration.....	33
6.1.1	Data Acquisition.....	35
6.2	System Shutdown.....	35

<b>Chapter 7:</b>	<b>Exposure .....</b>	<b>37</b>
7.1	Exposure Time .....	37
7.2	CCD Temperature .....	38
7.3	Dark Charge .....	39
7.3.1	Clean Until Trigger .....	41
7.3.1.1	Normal Shutter Mode .....	42
7.3.1.2	Open Before Trigger Mode .....	42
7.4	Saturation .....	43
<b>Chapter 8:</b>	<b>Analog to Digital Conversion .....</b>	<b>45</b>
8.1	Digitization .....	45
8.1.1	Speed .....	46
8.1.2	ADC Offset (Bias) .....	46
8.1.2.1	Correct Pixel Bias .....	46
8.2	Analog Gain .....	47
<b>Chapter 9:</b>	<b>Full Frame Readout .....</b>	<b>49</b>
9.1	Calculating Image Acquisition/Readout Time .....	51
9.1.1	CCD Readout Time .....	51
9.2	Experiment Timing .....	51
9.2.1	Trigger Response .....	53
9.2.1.1	No Response .....	53
9.2.1.2	Start on Single Trigger .....	56
9.2.1.3	Readout Per Trigger .....	58
9.2.1.4	Expose During Trigger Pulse .....	62
9.2.2	Trigger Determined By .....	65
9.3	Trigger Out .....	66
<b>Chapter 10:</b>	<b>Binning .....</b>	<b>67</b>
10.1	Hardware Binning .....	67
10.2	Software Binning .....	69
<b>Chapter 11:</b>	<b>Shutter Control .....</b>	<b>71</b>
11.1	Configuration .....	71
11.1.1	Mode .....	71
11.1.2	Opening Delay .....	72
11.1.3	Closing Delay .....	72
11.2	Using an External Shutter .....	73
<b>Appendix A:</b>	<b>Technical Specifications .....</b>	<b>75</b>
A.1	System Dimensions and Weight .....	75
A.1.1	Vacuum Window .....	75
A.2	Camera Specifications .....	75
A.2.1	Thermal Characteristics .....	76
A.3	Power Specifications .....	77
A.4	Environmental Specifications .....	78
A.4.1	Ventilation .....	78
A.5	External Chiller Specifications .....	78
A.6	Minimum Host Computer Specifications .....	78
<b>Appendix B:</b>	<b>Outline Drawings .....</b>	<b>81</b>
<b>Appendix C:</b>	<b>Drain Coolant from Lansis .....</b>	<b>85</b>
<b>Appendix D:</b>	<b>Custom Modes .....</b>	<b>87</b>

D.1	Custom Sensor .....	87
D.1.1	Custom Timing .....	88
<b>Appendix E: Troubleshooting .....</b>		<b>89</b>
E.1	General Camera Faults/Errors .....	89
E.1.1	Connection Failure or Logic Power Supply Overcurrent. ....	89
E.1.2	Power Supply Switch in On Position, But Power LED Extinguished ..	89
E.1.3	Overexposed CCD. ....	90
E.1.4	Camera Stops Working. ....	90
E.2	Cooling Faults/Errors. ....	91
E.2.1	Temperature Lock Cannot be Achieved or Maintained .....	91
E.2.2	Camera Loses Temperature Lock. ....	92
E.2.3	Gradual Deterioration of Cooling Capability .....	92
E.2.4	External Chiller: Low Coolant (Air in the Hoses) .....	92
E.3	LightField Faults/Errors .....	93
E.3.1	Devices Missing .....	93
E.3.2	Device is Occupied .....	94
E.3.3	Acquisition Started but Viewer Display Does Not Update .....	94
<b>Warranty and Service .....</b>		<b>97</b>
Limited Warranty .....		97
Basic Limited One (1) Year Warranty .....		97
Limited One (1) Year Warranty on Refurbished or Discontinued Products ..		97
XP Vacuum Chamber Limited Lifetime Warranty. ....		97
Sealed Chamber Integrity Limited 12 Month Warranty. ....		98
Vacuum Integrity Limited 12 Month Warranty .....		98
Image Intensifier Detector Limited One Year Warranty .....		98
X-Ray Detector Limited One Year Warranty .....		98
Software Limited Warranty. ....		98
Owner's Manual and Troubleshooting. ....		99
Your Responsibility .....		99
Contact Information .....		100



# Chapter 1: About this Document

---

Thank you for purchasing a Lansis camera system from Teledyne Princeton Instruments. Since 1981, Teledyne Princeton Instruments has been the legendary name behind the most revolutionary spectroscopy and imaging products for cutting edge research.

Please read this manual carefully before operating the camera. This will help you optimize the many features of this camera to suit your research needs.

If you have any questions about the information contained in this manual, contact the Teledyne Princeton Instruments customer service department. Refer to [Contact Information](#) on page 100 for complete contact information.

## 1.1 Intended Audience

This user manual is intended to be used by scientists and other personnel responsible for the installation, setup, configuration, and acquisition of imaging data collected using a Lansis system.

This document provides all information necessary to safely install, configure, and operate the Lansis, beginning with the system's initial installation.

## 1.2 Related Documentation

[Table 1-1](#) provides a list of related documentation and user manuals that may be useful when working with the Lansis camera system. To guarantee up-to-date information, always refer to the current release of each document listed.

**Table 1-1: Related Documentation**

Document Number	Document Title
–	LightField 6.15 (and newer) Online Help
–	Lansis Camera System Data Sheet
Varies	Spectrograph User Manual
Varies	PI Cam Programmer's Manual

## 1.3 Document Organization

This manual includes the following chapters and appendices:

- [Chapter 1, About this Document](#)  
This chapter provides information about the organization of this document, as well as related documents, safety information, and conventions used throughout the manual.
- [Chapter 2, Lansis Camera System](#)  
This chapter provides information about the components included with a standard Lansis camera system, as well as options that are available for purchase from Teledyne Princeton Instruments.
- [Chapter 3, Install LightField](#)  
This chapter provides information about the installation of Teledyne Princeton Instruments' LightField image acquisition software.
- [Chapter 4, System Block Diagrams](#)  
This chapter provides information about integrating the Lansis into a user's experiment.
- [Chapter 5, Hardware Configuration](#)  
This chapter provides information about the installation and configuration of system hardware.
- [Chapter 6, LightField First Light](#)  
This chapter provides a step-by-step procedure for placing a Lansis camera system in operation for the first time when using Teledyne Princeton Instruments' LightField 64-bit data acquisition software.
- [Chapter 7, Exposure](#)  
This chapter discusses the various factors that affect the signal acquired on the array, including array architecture, exposure time, temperature, and saturation.
- [Chapter 8, Analog to Digital Conversion](#)  
This chapter discusses the configuration of the Analog to Digital Conversion configuration parameters.
- [Chapter 9, Full Frame Readout](#)  
This chapter discusses Full Frame Readout operation and related parameter configuration.
- [Chapter 10, Binning](#)  
This chapter discusses the configuration of hardware and software binning.
- [Chapter 11, Shutter Control](#)  
This chapter discusses the configuration of shutter control parameters.
- [Appendix A, Technical Specifications](#)  
Provides CCD, system, and other basic specifications for a Lansis system.
- [Appendix B, Outline Drawings](#)  
Provides outline drawings of the Lansis camera and power supply.



- [Appendix C, Drain Coolant from Lansis](#)  
This appendix provides information necessary to safely drain coolant from within the Lansis camera body.
- [Appendix D, Custom Modes](#)  
Provides information necessary to configure custom chip modes on the Lansis.
- [Appendix E, Troubleshooting](#)  
This appendix provides recommended troubleshooting information for issues which may be encountered while working with a Lansis camera system.
- [Warranty and Service](#)  
This section provides warranty information for the Lansis. Contact information is also provided.

## 1.4 Safety Related Symbols Used in this Manual

The following safety symbols are used throughout this manual.



### **CAUTION!**

A **Caution** provides detailed information about actions and/or hazards that may result in damage to the equipment being used, including but not limited to the possible loss of data.



### **WARNING!**

A **Warning** provides detailed information about actions and/or hazards that may result in personal injury or death to individuals operating the equipment.



### **WARNING! RISK OF ELECTRIC SHOCK!**

The use of this symbol on equipment indicates that one or more nearby items pose an electric shock hazard and should be regarded as potentially dangerous. This same symbol appears in the manual adjacent to the text that discusses the hardware item(s) in question.

## 1.5 Lansis Safety Information

Before turning on the power supply, the ground prong of the power cord plug must be properly connected to the ground connector of the wall outlet. The wall outlet must have a third prong, or must be properly connected to an adapter that complies with these safety requirements.



### **WARNINGS!**

1. If the Lansis camera system is used in a manner not specified by Teledyne Princeton Instruments, the protection provided by the equipment may be impaired.
2. If the equipment or the wall outlet is damaged, the protective grounding could be disconnected. Do not use damaged equipment until its safety has been verified by authorized personnel. Disconnecting the protective earth terminal, inside or outside the apparatus, or any tampering with its operation is also prohibited.

Inspect the supplied power cord. If it is not compatible with the power socket, replace the cord with one that has suitable connectors on both ends.



### **WARNING!**

Replacement power cords or power plugs must have the same polarity and power rating as that of the original ones to avoid hazard due to electrical shock.



### **WARNING!**

If the Lansis is used in an OEM configuration in which the power supply is furnished by the OEM system, it is the responsibility of the OEM to ensure safety in accordance with good engineering practices.

## 1.6 Precautions

To prevent permanently damaging the Lansis system, observe the following precautions at all times.



### **CAUTION!**

---


1. The CCD array is very sensitive to static electricity. Touching the CCD can destroy it. Operations requiring contact with the device can only be performed at the factory.
  2. When using high-voltage equipment (e.g., an arc lamp,) with the camera system, be sure to turn the camera power ON LAST and turn the camera power OFF FIRST.
  3. Use caution when triggering high-current switching devices near the system (e.g., an arc lamp.) The CCD can be permanently damaged by transient voltage spikes. If electrically noisy devices are present, an isolated, conditioned power line or dedicated isolation transformer is highly recommended.
  4. Do not block air vents on the camera. Preventing the free flow of air overheats the camera and may damage it.
  5. If the Lansis camera system is used in a manner not specified by Teledyne Princeton Instruments, the protection provided by the equipment may be impaired.
-

*This page is intentionally blank.*

# Chapter 2: Lansis Camera System

This chapter provides an introduction to, and overview information about, Princeton Instrument's Lansis camera system. [Figure 2-1](#) shows those items that are typically included as part of a standard Lansis Camera system.

**Figure 2-1: Typical Lansis System Components**



4411-0174\_0001

Standard items for a typical air-cooled system include:

- Lansis Camera;
- Power Supply and Cable;
- USB3 Interface Cable;
- MCX to BNC Adapter Cables<sup>a</sup>;
- Certificate of Performance;
- Data Acquisition Software, and Hardware Key

a. Length May Vary

Standard items for a typical liquid-cooled system include:

- Lansis Camera;
- Power Supply and Cable;
- USB3 Interface Cable;
- MCX to BNC Adapter Cables<sup>a</sup>;
- External Chiller<sup>b</sup>;
- Coolant Hoses<sup>b</sup>;
- Certificate of Performance;
- Data Acquisition Software and Hardware Key.

a. Length May Vary

b. Not illustrated in [Figure 2-1](#).

## 2.1 Lansis Camera

Lansis cameras, illustrated in [Figure 2-2](#), are fully integrated camera systems. The camera contains all of the electronics necessary to read out and control the CCD device. For instance, it houses precision analog-to-digital converters (ADCs) positioned close to the CCD for lowest noise and has USB 3.0 electronics to interface with the host computer.

**Figure 2-2: Typical Lansis Cameras**



4411-0174\_0002

Lansis camera systems offer all basic CCD camera functions, such as Region of Interest (ROI) selection and binning, all under software control. It also provides advanced triggered operation as well as programmable TTL output.

Among the many state-of-the-art features are its maintenance-free permanent vacuum, integrated controller, deep thermoelectric air-cooling, and compact design.

### 2.1.1 eXcelon®

eXcelon is a CCD/EMCCD sensor technology jointly developed by Teledyne Princeton Instruments, Teledyne e2v, and Teledyne Photometrics®. Spectroscopy CCDs using this technology provide the following significant benefits:

- Improved Sensitivity  
Improved QE over broader wavelength region compared with back-illuminated sensors.
- Reduced Etaloning  
Up to 10 times lower etaloning or unwanted fringes in near infrared (NIR) region compared with standard back-illuminated CCDs.

## 2.1.2 Power

All voltages required by Lansis camera systems are generated and delivered by an external power supply included with each Lansis camera.



### CAUTION!

Use of a power supply other than that provided with the Lansis camera will void the camera warranty. For specific power supply requirements, contact Teledyne Princeton Instruments. Refer to [Contact Information](#) on page 100 for complete information.

The receptacle on the power supply should be compatible with the line-voltage line cords in common use in the region to which the system is shipped. If the power supply receptacle is incompatible, a compatible adapter should be installed on the line cord, taking care to maintain the proper polarity to protect the equipment and assure user safety.

[Figure 2-3](#) shows the connectors and indicators found on the rear of the Lansis power supply.



### REFERENCES:

Refer to [Section A.3, Power Specifications](#), on page 77 for detailed voltage specifications.

Figure 2-3: Lansis Camera Connectors and Indicators



4411-0174\_0079

Refer to [Table 2-1](#) for information about each connector and indicator on the power supply.

### 2.1.3 CCD Arrays

The Lansis camera system incorporates a back-illuminated, scientific-grade CCD to ensure the highest image fidelity, resolution, and acquisition flexibility required for scientific imaging. Large full wells, square pixels, and 100% fill factors provide high dynamic range and excellent spatial resolution. Your choice of CCD is already installed in the camera that you received and has been individually tested.

For complete specifications and information about the CCD used in Lansis cameras, refer to [Table A-3, Lansis CCD Specifications](#), on page 75.

### 2.1.4 Cooling

Dark current is reduced in Lansis camera systems by cooling the CCD array using Teledyne Princeton Instruments' exclusive ArcTec™ technology using air and/or circulating coolant. To prevent condensation and contamination from occurring, cameras cooled this way are evacuated. Refer to [Table A-4, Default Operating Temperature](#), on page 76 for specific cooling information.

#### 2.1.4.1 Internal Fan

The Lansis camera is equipped with an internal cooling fan that:

- Removes heat from the Peltier device that cools the CCD array, and
- Cools the electronics.

Teledyne Princeton Instruments ArcTec cooling technology cools Lansis's CCD assisted by air drawn into the camera by an internal fan mounted on the rear of the camera. The circulating air then vents out through slots on the side panels. By default, the fan is always in operation and air-cooling of both the CCD and the internal electronics occurs continuously. In most cases, the low-vibration fan action does not adversely affect image acquisition. However, in some applications, the fan's vibration could reduce image quality. In these instances, the internal fan can be disabled on the Sensor Expander within LightField.

- When Lansis is being air cooled (i.e., no external chiller is in use,) the internal fan can be disabled for a brief period (i.e., a few seconds.) Sensors within Lansis monitor its temperature and will reactivate the fan before temperatures within the camera rise too much.
- When Lansis is being cooled using an external chiller, the internal fan can be safely disabled for a longer period. As with air cooled applications, internal sensors monitor the temperature and will reactivate the fan if necessary.

For the fan to function properly, uninhibited air circulation must be maintained between the sides of the camera and the laboratory atmosphere.

#### 2.1.4.2 External Cooling Circulator

Lansis cameras equipped for liquid cooling can provide a low vibration system for data acquisition when connected to a chiller. Any commercially available circulator can be used provided it is capable of continuously pumping a 50:50 mixture of room temperature (23°C) water and ethylene glycol antifreeze at 1 liter per minute. Refer to [Section A.5, External Chiller Specifications](#), on page 78 for additional information. If desired, contact Teledyne Princeton Instruments for additional recommendations. Refer to [Contact Information](#) on page 100 for complete information.



**CAUTION!**

Never set the coolant temperature below the dew point.

Coolant temperature should not be below the dew point of the ambient air. Internal condensation caused by operation below the dew point may damage the camera and will void the warranty. Additionally, the Lansis monitors its internal temperature and will automatically restart its fan regardless of the software settings if it senses excessive internal heat buildup.

**2.1.4.3 Coolant Ports**

The liquid cooled model of the Lansis camera is equipped with cooling ports that allow it to be connected to an external chiller. As is the case with circulating air, circulating coolant removes heat produced by the camera. This means of heat removal is designed for vibration-free data acquisition. For the circulating coolant to function properly, free air circulation must be maintained between the sides of the external circulator and the laboratory atmosphere.

**CAUTION!**

NEVER apply negative pressure to the liquid circulator fittings on the Lansis camera. Doing so may permanently damage the camera.

**CAUTION!**

The wetted areas of the Lansis heat exchanger and fittings are nickel plated.

2.1.5 Rear-Panel Connectors and Indicators

Figure 2-4 illustrates the rear-panel connectors and indicators on a Lansis camera.

Figure 2-4: Lansis Rear-Panel Connectors and Indicators



4411-0174\_0003

Refer to Table 2-1 for information about each rear-panel connector and indicator.

Table 2-1: Lansis Rear-Panel Connectors and Indicators

Label	Description
USB3	Control signals and data are transmitted between the camera and host computer via this port.
TRIGGER IN	0 – +5.0 V <sub>DC</sub> TTL-compatible logic level input with a 10 kΩ pull-up resistor. Allows data acquisition and readout to be synchronized with external events. Positive or negative edge triggering is programmable. <b>CAUTION:</b> Do not apply a voltage greater than +5.0 V <sub>DC</sub> or less than 0 V <sub>DC</sub> (i.e., a negative voltage,) to this input. Doing so may permanently damage the Lansis camera.
OUT1	0 to +5.0 V <sub>DC</sub> programmable TTL-compatible logic level output. This output can be programmed and inverted via the application software
OUT2	0 to +5.0 V <sub>DC</sub> programmable TTL-compatible logic level output. This output can be programmed and inverted via the application software.
POWER	Power input from external power supply provided with the Lansis system.
STATUS	This LED indicates the state of the camera. See Table 2-2 to decode the camera state.



**Table 2-2: Lansis Rear Panel LED Indications**

LED	Lansis system state
Green, solid	No errors, camera is in standby (executing cleans)
Off	Acquiring-No errors, camera is in acquisition.
Yellow, Blinking	The TEC heat sink is too hot, (hot-side over-temperature error), the FPGA turned off the TEC Supply. This should allow the heat sink to cool. If/when the heatsink cools enough, the FPGA will set the LED back to green or off and re-enable the TEC supply.
Red, Blinking	Critical Board Temperature (TEC remains driven). This will automatically clear if the board temperature goes down to an acceptable level.
Yellow, Solid	TEC power supply has a failure.
Red, Solid	Other power supply has a failure.

## 2.2 Cables

[Table 2-3](#) describes the cables included with a standard Lansis Camera System.

**Table 2-3: Standard Lansis Camera System Cables**

Cable	Part Number	Description/Purpose	Length
USB 3	6050-0738	Connects the USB 3 connector on the rear of the Lansis with the host computer. 	~5 m
MCX to BNC	6050-0540	Three (3) MCX to BNC adapter cables are included. These allow connection to the MCX connectors on the rear of the Lansis using BNC cables.. 	~15cm

## 2.3 Certificate of Performance

Each Lansis camera is shipped with a Certificate of Performance which states that the camera system has been assembled and tested according to approved Teledyne

**EAR99 Technology Subject to Restrictions Contained on the Cover Page.**

Princeton Instruments procedures. It documents the camera's performance data as measured during the testing of the Lansis and lists the following camera- and customer-specific information:

- Sales Order Number;
- Purchase Order Number;
- Camera Serial Numbers

This information is useful when contacting Teledyne Princeton Instruments Customer Support.

## 2.4 Application Software

Teledyne Princeton Instruments offers a number of data acquisition software packages for use with Lansis camera systems, including:

- **LightField®**  
The Lansis camera can be operated using LightField, Teledyne Princeton Instruments' 64-bit Windows® compatible software package. LightField combines complete control over Teledyne Princeton Instruments' cameras and spectrographs with easy-to-use tools for experimental setup, data acquisition and post-processing. LightField makes data integrity priority #1 via automatic saving to disk, time stamping and retention of both raw and corrected data with full experimental details saved in each file. LightField works seamlessly in multi-user facilities, remembering each user's hardware and software configurations and tailoring options and features accordingly. The IntelliCal® package is the highest-performance wavelength calibration software available, providing up to 10X greater accuracy across the entire focal plane than competing routines.
- **PICam™**  
The standard 64-bit software interface for cooled CCD cameras from Teledyne Princeton Instruments. PICam is an ANSI C library of camera control and data acquisition functions. PICam is available for 64-bit Windows and Linux systems. Refer to the PICam Programmer's Manual for the list of supported operating systems.

## 2.5 Accessories

Teledyne Princeton Instruments offers a number of optional accessories that are compatible with Lansis. For complete ordering information, contact Teledyne Princeton Instruments.

### 2.5.1 External Shutters

The LANSIS camera does not directly drive electro-mechanical shutters. A logic signal indicative of the time the shutter should be open is available at the logic output connectors if it is selected in the software. This signal can be used to activate an external shutter driver. Contact Teledyne Princeton Instruments for cases in which use of an external shutter is necessary.

## 2.6 Unpack the System

All required items should be included with the shipment. The Lansis system has been manufactured according to the camera options specified at the time of purchase, including the CCD window and coatings that were ordered.

When unpacking the system, examine the system components for any signs of shipping damage. If there are any, notify Teledyne Princeton Instruments immediately and file a claim with the carrier. Be sure to save the shipping carton for inspection by the carrier. If damage is not apparent but system specifications cannot be achieved, internal damage may have occurred in shipment.

Retain all original packing materials so that the Lansis system can be easily and safely packaged and shipped to another location or returned for service if necessary. If assistance is required at any time, contact Teledyne Princeton Instruments Customer Support. Refer to [Contact Information](#) on page 100 for complete information.

### 2.6.1 Verify Equipment and Parts Inventory

Verify all equipment and parts required to set up the Lansis system have been delivered.

A typical system consists of:

- Lansis Camera;
- Power Supply and Cable;
- USB3 Interface Cable;
- MCX to BNC Adapter Cables<sup>a</sup>;
- Certificate of Performance;
- Data Acquisition Software and Hardware Key

a. Length May Vary

Accessories that may have been purchased include:

- External Chiller and hoses (for liquid-cooled system).

## 2.7 Lansis Camera and System Maintenance



### **WARNING!**

Turn off all power to the equipment and secure all covers before cleaning the units. Otherwise, damage to the equipment or injury to you could occur.

### 2.7.1 Camera

Although there is no periodic maintenance that needs to be performed on a Lansis camera, users are advised to wipe it down with a clean damp cloth from time to time. This operation should only be done on the external surfaces and with all covers secured. In dampening the cloth, use clean water only. No soap, solvents or abrasives should be used. Not only are they not required, but they could damage the finish of the surfaces on which they are used.

### 2.7.2 Optical Surfaces

As a good practice, the camera must be closed/capped off with the supplied dust cover or lens cap when not in use. If it becomes necessary to clean the optical window due to the accumulation of atmospheric dust, we advise that the drag-wipe technique be used. This involves dipping a clean cellulose lens tissue into clean anhydrous methanol, and then dragging the dampened tissue over the optical surface to be cleaned. Do not allow any other material to touch the optical surfaces.

### 2.7.3 Repairs

Because the Lansis camera system contains no user-serviceable parts, repairs must be performed by Teledyne Princeton Instruments. Should the system need repair, contact Teledyne Princeton Instruments customer support for instructions. Refer to [Contact Information](#) on page 100 for complete information.

Save the original packing materials and use them whenever shipping the system or system components.

# Chapter 3: Install LightField

---

This chapter provides the installation procedure for LightField application software.



---

**NOTE:**

If LightField has already been successfully installed on the host computer, this chapter may be skipped.

---

## 3.1 Prerequisites

Before beginning to install LightField, verify that:

- The operating system on the desired host computer is Windows® 7/8/10 (64-bit);
- The host computer supports USB3;  
If it does not support USB3, refer to the host computer manufacturer's instructions for installing a USB3 interface card;
- The installation media and hardware key are available.

## 3.2 Installation Procedure

Perform the following procedure to install LightField on the host computer:

1. Insert the LightField Installation media on the host computer, and follow the on-screen prompts.
2. After the installation has been completed, reboot the host computer.
3. Connect the Lansis system components to the host computer and apply power.
4. Launch LightField, activate it, and begin experiment configuration.


*This page is intentionally blank.*



# Chapter 4: System Block Diagrams


This section provides block diagrams of typical system configurations.

Figure 4-1: Block Diagram: Typical Air-Cooled Experiment



4411-0174\_0013

Figure 4-2: Block Diagram: Typical Liquid-Cooled Experiment



*This page is intentionally blank.*

# Chapter 5: Hardware Configuration

---

This chapter provides information about the installation of hardware that may be installed on a Lansis camera.

A Lansis camera must be properly mounted to a spectrograph in order to focus the optics. Additional precautions must also be taken to prevent overexposure of the camera.

The distance to the focal plane from the front of the mechanical assembly depends on the specific configuration.

## 5.1 Spectrograph Adapters

Teledyne Princeton Instruments offers a number of adapters for mounting a Lansis camera onto a spectrograph, including:

- Teledyne Acton Research Series Spectrographs, and
- IsoPlane Family of Spectrographs.



---

**NOTE:**

Other adapters may be available. Contact Teledyne Princeton Instruments for specific information. Refer to [Contact Information](#) on page 100, for complete contact information.

---

## 5.2 Overexposure Protection

Cameras that are exposed to room light or other strong light sources can become saturated. This most often occurs when operating without a shutter or with the shutter locked open. When LANSIS is mounted to a spectrograph with high incoming light, close the entrance slit of the spectrograph to a normal spectroscopy setting to reduce the incident light.



---

**NOTE:**

If the CCD is cooled to low temperatures (below -50°C), exposure to ambient light can over-saturate it, causing charge to be stored in the bulk silicon. This charge is released over time, and manifests as increased dark current. If the camera shows elevated dark charge, to restore dark charge to its original level it may be necessary to bring the camera back to room temperature, then cool it down again in the dark. At low temperatures, it takes a long time for the CCD to return to its normal dark current level. The CCD is not permanently harmed by over-saturation

---

## 5.3 External Chiller Use

For liquid-cooled cameras, an external chiller provides a vibration-free method of heat removal. Perform the following procedure to connect an external chiller to a liquid-cooled Lansis camera:

**NOTE:**

For specific configuration information, refer to the manufacturer-supplied documentation included with the external chiller.

1. Verify the camera and the circulator power switches are turned off.
2. Verify the chiller is positioned a minimum of 6 inches (150 mm) below the Lansis camera. The vertical distance should not exceed 10 feet (3 m).  
Typically, the camera is at table height and the circulator is on the floor.
3. Using the supplied cooling hoses, connect the two coolant ports on the external chiller to the two coolant ports on Lansis.

**CAUTION!**

Verify there are no kinks in hoses that may impede coolant flow. Lack of sufficient flow can seriously harm the detector. Any resulting damage is not covered under warranty.

**NOTE:**

Although there are no dedicated IN or OUT coolant ports on a Lansis camera, Teledyne Princeton Instruments recommends that the lower port (depending on the camera's physical orientation,) be used for the coolant inlet for best efficiency.

**NOTE:**

Damage caused by water leaking into the Lansis voids the warranty.

4. Verify that the reservoir on the external chiller contains sufficient coolant as specified by its manufacturer. If additional coolant is required, use a 50:50 mixture of water and ethylene glycol based commercial antifreeze to add sufficient coolant.
5. Replace the reservoir cap.
6. Plug the external circulator into a 100-240 V<sub>AC</sub>, 47-63 Hz power source.
7. Turn the chiller on.

**CAUTION!**

When configuring an external chiller, adhere to the following:

- Coolant flow-rate should be at least 3 l/min (0.8 gpm)
- Coolant pressure should never exceed 137kPa (20 psi).
- Never set the circulator temperature below the dew point.

8. Verify there are no leaks or air bubbles in the hoses.

**NOTE:**

Small air bubbles (approximately the size of bubbles in soda) are common, particularly immediately following start up. These bubbles do not prevent proper operation.

- If no problems are observed, proceed to step 9.
  - If there are leaks or air bubbles, turn the circulator off and correct the problem(s) by securing the hoses or adding more coolant to the reservoir. Turn the circulator back on. Recheck and if there are no problems, proceed to step 9.
9. Turn on the camera.
  10. Launch LightField.

### 5.3.1 Experiment Shutdown

Following the completion of a liquid-cooled experiment, perform the following procedure:

1. Turn off the Lansis camera's power supply.
2. Turn off the chiller according to all manufacturer-supplied documentation and procedures.
3. If desired, carefully disconnect the coolant hoses from the Lansis camera.

**WARNING!**

If the Lansis is to be shipped, to avoid potential catastrophic damage to the camera, all coolant must be drained from it. Refer to [Appendix C, Drain Coolant from Lansis](#), on page 85.

*This page is intentionally blank.*

# Chapter 6: LightField First Light


Once the Lansis camera has been configured as described in [Chapter 5, Hardware Configuration](#), acquiring data using LightField is straightforward. For most applications, simply:

- Establish optimum performance using Preview mode;
- Set a target camera temperature;
- Wait until the system's temperature has stabilized;
- Acquire live data in Acquire mode.

Additional considerations regarding experiment setup and equipment configuration are addressed in the LightField Online Help.

During data acquisition, the CCD array is exposed to a source and charge accumulates in the pixels. After the defined exposure time, the accumulated signal is readout of the array, digitized, and then transferred to the host computer. Upon data transfer, the data are displayed and/or stored via the application software. This sequence is illustrated by the block diagram shown in [Figure 6-1](#).

**Figure 6-1: Block Diagram for Lansis Systems**



4411-0174\_0020

Whether or not the data are displayed and/or stored depends on the data collection operation that has been selected in the application software:

- **Preview**

Data collection operations use the Experiment Setup parameters to establish the exposure time (the period when signal of interest is allowed to accumulate on the CCD). Preview is typically used when setting up the system

In Preview mode, the number of frames is ignored. A single frame is acquired and displayed, another frame is acquired which overwrites the currently displayed data, and so on, until Stop is selected.

**NOTE:**

The last frame acquired before Stop is selected cannot be stored.

However, the last frame of data can be saved by configuring the Number of Frames to 1 and clicking Acquire rather than Stop.

Preview mode is particularly convenient for familiarization and configuration. For ease in focusing, the screen refresh rate should be as quick as possible, achieved by operating with axes and cross-sections off, and with Zoom 1:1 selected.

- **Acquire**

Acquire mode is typically used for the collection and storage of data. In Acquire mode, every frame of data collected can be automatically stored, so the completed dataset may include multiple frames with one or more accumulations. This mode is typically selected during actual data collection.

One limitation of Acquire mode operation is that if data acquisition continues at too fast a rate for it to be stored, data overflow may occur. If this happens, try the following to remedy the situation:

- Shut down any other programs that may be running and utilizing resources;
- Decrease the camera speed;
- Use a faster host computer.



## 6.1 Set Up and Configuration

This section provides step-by-step instructions for acquiring a spectrum in LightField for the first time. The intent of this procedure is to gain familiarity with the operation of the system and to show that it is functioning properly. Once basic familiarity has been established, additional, more complex configurations can be implemented.

The following procedure assumes:

- The system has been set up in accordance with the instructions in previous chapters;
- Familiarity with LightField;  
If this is not the case, refer to the online help while performing this procedure.
- A suitable light source (e.g., mercury pen-ray lamp,) has been mounted in front of the entrance slit of the spectrograph.

Any light source with line output can be used. Standard fluorescent overhead lamps have good calibration lines as well. If there are no line sources available, it is possible to use a broadband source such as tungsten for the alignment. If this is the case, use a wavelength setting of 0.0 nm for alignment purposes.




### CAUTION!

**Overexposure Protection:** Cameras that are exposed to room light or other continuous light sources will quickly become saturated. If the camera is mounted to a spectrograph, close the entrance slit of the spectrograph to 20 mm (typical) to reduce the incident light.

Perform the following procedure to set up and configure the system to acquire a spectrum:

1. Set the spectrograph entrance slit width to minimum (20 mm if possible.)
2. Turn ON the spectrograph (if applicable.)
3. Mount a light source such as a Teledyne Princeton Instruments IntelliCal<sup>®</sup> Hg/Ne-Ar Dual Switchable light source in front of the entrance slit.
4. Power on the camera.
5. Launch LightField.
6. Once LightField has started, a Lansis camera icon, will be shown in the Available Devices area. See [Figure 6-2](#).

Figure 6-2: Typical Available Devices Area




4411-0174\_0051

7. Drag the camera icon into the Experiment Devices area.

The Experiment Settings stack on the left includes several expanders. Since this is a new experiment, the default configuration settings for the camera are used. See [Figure 6-3](#).

**Figure 6-3: Experiment Devices Area**



The Status bar across the bottom of the window includes an icon for temperature status. Temperature Status reports the current system temperature and whether the set temperature has been reached. Clicking on the icon opens the Sensor expander in which the desired temperature can be configured.



#### **NOTE:**

Each Lansis camera is carefully calibrated at the factory for optimum performance at low CCD temperatures. If operated at warmer temperatures, such as during system cool down or if the temperature set point is set warmer than the factory-calibrated temperature, the camera might not image properly. Artifacts in the image, fixed value pixels, or lack of imaging in one or more regions may occur.

8. Set the spectrometer to a grating and wavelength that will produce a visible spectrum with the light source in use.



#### **NOTE:**

Overhead fluorescent lights produce a mercury spectrum. Use a white card tilted at 45 degrees in front of the entrance slit to reflect overhead light into the spectrograph. Select 500 nm as the spectral line.

9. Turn on the light source at the spectrograph's entrance slit.
10. Click Run to begin previewing the data.  
Depending on the display settings, you should see either a spectral band (image) or a graph. Background noise will decrease as the camera cools to its default temperature.
11. Turn off the light source. The data display should change to a background noise pattern or low intensity graph.  
If this occurs you have confirmed that light entering the spectrograph is being seen by the camera.
12. If there is little or no difference between the data displayed when the light source is on or off:
  - a. Verify that the light source has power and is turned on.
  - b. Verify that the entrance slit is open at least 10 mm.
  - c. Check the Exposure Time (Common Acquisition Settings expander).
  - d. Verify that there is no shutter between your spectrograph and the camera. If there is a shutter and you wish to drive it using TTL output signals from the Lansis, verify that the cables are connected properly and the camera output signals are correctly configured.
13. For additional assistance, contact Customer Support. Refer to [Contact Information](#) on page 100 for complete information.

### 6.1.1 Data Acquisition

Perform the following procedure to acquire live data:

1. After the system has been aligned and focused, stop the camera.
2. Make any required changes to the experiment setup and software parameters.  
Changes may include:
  - Adjusting the exposure time;
  - Setting up an entrance slit shutter;
  - Changing the timing mode to External Sync;
  - Lowering the temperature.
3. Begin running Acquire mode. Data will be acquired and displayed/stored based on the experiment settings.
4. Once data acquisition is complete, either:
  - Leave the camera power on so the array temperature will remain locked for subsequent data acquisition;
  - Shut down the system. Refer to [Section 6.2, System Shutdown](#), for proper shutdown procedures.

## 6.2 System Shutdown

Perform the following procedure to shutdown the Lansis system:

1. Set the camera's temperature to 0 °C and allow it to come to temperature.
2. Exit LightField.
3. Turn off the light source.
4. Turn off the spectrograph.
5. Turn off the camera power.

*This page is intentionally blank.*

# Chapter 7: Exposure

This chapter discusses factors that may affect the signal acquired on the CCD array, such as:

- Exposure Time;
- CCD Temperature;
- Dark Charge;
- Saturation.

## 7.1 Exposure Time

Exposure time is the time between commands sent by the data acquisition software to start and stop signal accumulation on the sensor.



### NOTE:

Exposure time is configured on the Common Acquisition Settings expander.


In combination with triggers these commands control when continuous cleaning of the CCD stops and when the accumulated signal will be read out.

Cleaning prevents the buildup of dark current and unwanted signal before the start of the exposure time. At the end of the exposure time, the CCD is readout and, depending on the specific experiment configuration:

- The next exposure begins (i.e., when configured for Free Run mode);
- Cleaning starts and continues until an incoming trigger is received (i.e., when configured for Trigger mode.)

Lansis requires an external shutter in order to control exposure of its CCD. [Figure 7-1](#) illustrates how the exposure period is measured. The Reading Out signal at the OUT 1 and OUT 2 connectors on the back of the Lansis can be used to monitor the exposure and readout cycle ( $t_R$ ) which is also shown in [Figure 7-1](#).

**Figure 7-1: Timing Diagram: CCD Exposure with Shutter Compensation**



where:

- $t_o$  is the shutter opening compensation time;  
This value is automatically configured for Lansis.
- $t_c$  is the shutter closing compensation time;
- $t_{exp}$  is the exposure time;
- $t_R$  is the readout time.

Note that the READOUT signal is:

- LOW during Shutter Open/Close Compensation Times;
- LOW during Exposure Time;
- HIGH during Readout.

Since the shutter behaves like an iris, the opening and closing of the shutter will cause the center of the CCD to be exposed slightly longer than the edges. It is important to realize this physical limitation, particularly when using short exposures.



### CAUTION!

A Shutter can overheat when short, rapidly repeated exposures are used, or if the shutter is held open for an extended period of time.

## 7.2 CCD Temperature

Lowering the temperature of the CCD generally enhances the quality of an acquired signal. Temperature control is configured on the Sensor expander.

Once a Temperature Setpoint has been programmed on the Sensor Expander, the software controls the camera's cooling circuits to reach the programmed array temperature.

When the camera's temperature has cooled to within  $\pm 0.5^{\circ}\text{C}$  of the programmed setpoint, Locked is displayed to the right of the target temperature in LightField's status bar. At this point, Lansis's control loop will continue to reduce the thermal error, typically holding the temperature to within  $\pm 0.05^{\circ}\text{C}$  of the programmed setpoint.

The time required to achieve lock may vary considerably depending on factors such as the camera type, CCD array type, ambient temperature, etc. Ultimate system stability is achieved approximately 20 minutes after thermal loop locks. However, as long as Lansis's CCD has cooled sufficiently to image, focusing of the system can begin. Typically, quantitative data may be acquired as soon as thermal lock is relatively stable. If, however, an experiment requires exposure times in the 10 to 30 minutes range or more, it is recommended that Lansis be permitted to cool for a longer period of time before acquiring live data.

The deepest operating temperature for a system depends on the CCD array size and packaging. Refer to [Table A-4, Default Operating Temperature](#) on page 76 for typical deepest cooling temperatures.

## 7.3 Dark Charge

Dark charge (or dark current) is the thermally induced buildup of charge in the CCD over time. The statistical noise associated with this charge is known as dark charge noise. Dark charge values vary widely from one CCD array to another and are exponentially temperature dependent. In the case of cameras with MPP type arrays, the average dark charge is extremely small. However, the dark-charge distribution is such that a significant number of pixels may exhibit a much higher dark charge, limiting the maximum practical exposure. Dark charge effect is more pronounced in the case of cameras having a non-MPP array (such as deep-depletion devices.)

With the light into the camera completely blocked, the CCD will collect a dark charge pattern, dependent on the exposure time and CCD temperature. The longer the exposure time and the warmer the camera, the larger and less uniform this background will appear. Thus, to minimize dark-charge effects, you should operate with the lowest CCD temperature possible.



### CAUTION!

If a sudden change in the baseline signal is observed, there may be excessive humidity in the camera vacuum enclosure. Turn off the camera and contact Teledyne Princeton Instruments Customer Support. Refer to [Contact Information](#) on page 100 for complete contact information.

Do not be concerned about the DC level of this background. This is not noise but rather a subtractable bias pattern. By acquiring/saving a dark charge background image under conditions identical to those that will be used to acquire live data, this background image can then be subtracted from the acquired image, thus reducing dark-charge effects.

Although the dark charge pattern will be subtracted from the acquired image, both the dark charge pattern and the acquired image include system readout noise,  $N_R$ , which is mean square additive. Therefore, when acquiring a dark charge pattern, it is strongly recommended that multiple frames of dark data be acquired and averaged since the cumulative noise within the dark background is reduced by the square root of the number of frames that have been averaged. The cumulative readout noise, which is the sum of the acquired image readout noise plus the dark data readout noise, is calculated as follows:

$$N_{RT} = \sqrt{N_R^2 + \frac{N_R^2}{F}}$$

where:

- $N_{RT}$  = Total Readout Noise;
- $N_R$  = Readout Noise;
- $F$  = Number of Dark Pattern Frames Acquired.

For example, when acquiring one frame of dark data to generate the subtractable pattern, the dark pattern readout noise will be  $\sim 1.414\sqrt{N_R}$ , which is an increase of  $\sim 41\%$  over the baseline readout noise. This value is calculated as follows:

$$\begin{aligned}
 N_{RT} &= \sqrt{N_R^2 + \frac{N_R^2}{F}} \\
 &= \sqrt{N_R^2 + \frac{N_R^2}{1}} \\
 &= \sqrt{N_R^2 + \frac{N_R^2}{1}} \\
 &= \sqrt{2N_R^2} \\
 &= 1.414N_R
 \end{aligned}$$

Refer to [Table 7-1](#) for typical noise penalty figures as a percentage of readout noise,  $N_R$ .

**Table 7-1: Readout Noise Penalty vs Number of Acquired Frames**

Number of Acquired Frames (F)	Noise Penalty <sup>a</sup>
1	41%
4	12%
8	6%
16	3%

a. Expressed as a percentage of  $N_R$



### 7.3.1 Clean Until Trigger

When using an external trigger to initiate data readout, Lansis supports Clean Until Trigger (CUT,) an additional level of cleaning/removing accumulated dark charge that continues until the moment the External Sync pulse is received.




#### REFERENCES:

For information about the use and configuration of external triggers, refer to:

- [Section 9.2, Experiment Timing](#), on page 51, for Full Frame Mode;

Figure 7-2 illustrates a flowchart of this mode.

Figure 7-2: Flowchart: Clean Until Trigger




4411-0174\_0026

### 7.3.1.1 Normal Shutter Mode

When an incoming Trigger pulse has been received, cleaning of the array stops as soon as the current cleaning pattern has been completed and shifted. Because the incoming trigger is not synchronous with the cleaning cycle, there is an inherent jitter of up to one cleaning cycle in the system's response to an incoming trigger. After this, the shutter remains open for the programmed exposure time, and frame collection can begin.

Figure 7-3 illustrates the timing diagram for Clean Until Trigger, Normal shutter mode. In this figure, system jitter is shown as light gray dashed lines.

Figure 7-3: Timing Diagram: Clean Until Trigger, Normal Shutter Mode




4411-0174\_0027

### 7.3.1.2 Open Before Trigger Mode

With Open Before Trigger shutter mode, the shutter is opened at the beginning of Clean Until Trigger. As with Normal mode, once the incoming Trigger pulse has been received, cleaning of the array stops as soon as the current cleaning pattern has been completed and shifted. Because the incoming trigger is not synchronous with the cleaning cycle, there is an inherent jitter of up to one cleaning cycle in the system's response to an incoming trigger. After this, the shutter remains open for the programmed exposure time, and frame collection can begin.

Figure 7-4 illustrates the timing diagram for Clean Until Trigger, Open Before Trigger shutter mode. In this figure, system jitter is shown as light gray dashed lines.

Figure 7-4: Timing Diagram: Clean Until Trigger, Open Before Trigger



4411-0174\_0099

If the trigger arrives while a cleaning pattern is executing, the pattern will be completed, then cleaning will cease. This could result in vertical smear. The cleaning pattern is under user control in LightField, and the use of a cleaning pattern with only one vertical shift will reduce the vertical smear to at most one row.

As expected, the response latency is on the order of one vertical shift time, from 1-30  $\mu\text{s}$  depending on the array. This latency does not prevent the incoming signal from being detected, since photo generated electrons are still collected over the entire active area. However, if the signal arrival is coincident with the vertical shifting, image smearing of up to one pixel is possible. The amount of smearing is a function of the signal duration compared to the single vertical shift time.

## 7.4 Saturation

When signal levels in some part of the image are very high, charge generated in one pixel may exceed the well capacity of the pixel, spilling over into adjacent pixels in a process called blooming. In this case a shorter exposure is advisable, with signal averaging to enhance S/N (Signal-to-Noise ratio) accomplished through the software.

For signal levels low enough to be readout-noise limited, longer exposure times, and therefore longer signal accumulation in the CCD, will improve the S/N ratio approximately linearly with the length of exposure time. However, due to cosmic ray strikes, the maximum practical exposure time is in the range of 15 to 30 minutes. Be aware that some deep depletion CCDs have a higher dark charge rate, may be dark current limited, and may also be more sensitive to cosmic rays. High- $\rho$  ( $\rho$ ) CCDs are markedly more sensitive to cosmic rays.

*This page is intentionally blank.*

# Chapter 8: Analog to Digital Conversion

After the programmed exposure time has elapsed, accumulated charge stored in the CCD array must be:


- Read out;
- Converted to a digital format;
- Transferred to the application software where it can be displayed and/or stored.

Analog-to Digital conversion factors are configured within LightField on the Analog to Digital Conversion expander. The following parameters are able to be configured:

- Speed;
- Analog Gain.

See [Figure 8-1](#).

Figure 8-1: Typical Analog to Digital Conversion Expander



The following sections describe the impact each of these parameters has on acquired image data.

## 8.1 Digitization

After gain has been applied to the signal, the Analog-to-Digital Converter (ADC) converts analog information (continuous amplitudes) into digital data (quantified, discrete steps,) that can be read, displayed, and stored by the application software. The number of bits per pixel is, by design, set to 16.

Factors associated with digitization include:

- **Speed**  
This is software configurable within LightField, and specifies the rate at which data are digitized.
- **ADC Offset (Bias).**  
This is not a user-configurable value; it is set at the factory at the time of manufacture.

These factors are discussed in the following paragraphs.

### 8.1.1 Speed

Because the readout noise for a CCD array increases with readout rate, it may be necessary at times to trade off readout speed for high dynamic range.

Lansis supports the following digitization speeds:

- 100 kHz
- 2 MHz

### 8.1.2 ADC Offset (Bias)

With the camera's light path completely blocked, the CCD accumulates a dark charge pattern, dependent on the exposure time and camera temperature. The longer the exposure time and the warmer the camera, the larger this background will appear. A 500-600 ADU offset is included to avoid loss of signal when noise and/or drift might inadvertently result in a signal  $< 0$  ADU. This offset value ensures that all the true variation in the signal can really be seen and not lost below the A/D "0" value. Since the offset is added to the signal, these counts only minimally reduce the range of the signal from 65535 (16-bit A/D) to a value in the range of 500-600 counts lower.



#### **WARNING!**

If a sudden change in the baseline signal is observed, there may be excessive humidity in the camera vacuum enclosure. Turn off the camera and contact Teledyne Princeton Instruments Customer Support. Refer to [Contact Information](#) on page 100 for complete information.




#### **NOTES:**

1. It is important to note that the bias level is not noise. It is a fully subtractable readout pattern. Every device has been thoroughly tested to ensure its compliance with Teledyne Princeton Instruments' demanding specifications.
2. The ADC Offset is pre-set at the factory and is not user-changeable.

#### 8.1.2.1 Correct Pixel Bias

By default, Pixel Bias Correction is ENABLED in LightField. Correct Pixel Bias automatically corrects pixel bias drift that may be introduced.

Pixel bias correction operates by reading an artificially created row, which has essentially no charge, at the start of the frame and subtracting the row average from all following pixels in the frame. Therefore it adds one row readout time to the frame time. This is particularly useful when data are to be collected over a long time and maximum baseline stability is important.

This setting can be disabled by clicking  on the Analog to Digital Conversion expander, and deselecting the Correct Pixel Bias option.

## 8.2 Analog Gain

Controller gain, a configurable function of the preamplifier, changes the relationship between the number of electrons acquired on the CCD and the Analog-to-Digital Units (ADUs) generated. The level of gain is configured by the Analog Gain parameter on the Analog to Digital Conversion expander.

Supported Analog Gain options are:

- Low

This gain setting is best suited for binning applications, particularly when high-level signals are being digitized. Although this setting may be used with non-binned modes, the CCD single pixel well typically will not reach ADC saturation when applying Low gain.

- High.

Applications that consistently measure low-level signals should select this gain setting since it requires fewer electrons to generate each ADU. Additionally, this setting can reduce some sources of noise.

**NOTE:**

The Certificate of Performance included with Lansis specifies the measured gain values for various Analog Gain settings.

Actual electron counts required are also dependent upon the configured readout rate.

**Table 8-1: Typical Electron Counts vs. Analog Gain Setting**

Low	High
5 e <sup>-</sup> /count	0.7 e <sup>-</sup> /count

*This page is intentionally blank.*



# Chapter 9: Full Frame Readout

When operating in Full Frame mode, Lansis reads and processes a complete frame of data at a time. Every pixel of information is digitized individually.

Figure 9-1 illustrates a CCD array following exposure but prior to the beginning of readout.



 **NOTE:** Capital letters represent different amounts of accumulated charge, both desired signal and dark charge.


Figure 9-1: Full Frame Readout: Unshifted CCD Charge



4411-0174\_0029

Readout of the CCD begins with the simultaneous shifting of all pixels one row toward the Shift Register which is a single row of pixels along the edge of the CCD. The Shift Register is not sensitive to light and is only used to store charge during readout. See [Figure 9-2](#).

Figure 9-2: Full Frame Readout: One Row of Charge Shifted into Shift Register



4411-0174\_0063




**NOTE:**

Typically, each pixel within the Shift Register can store up to twice as much charge as pixels within the CCD imaging area.

Once one row of charge has been shifted into the Shift Register, the charge is then shifted by one pixel toward the corresponding Output Node where it is then digitized. See [Figure 9-3](#).

Figure 9-3: Full Frame Readout: One Pixel of Charge Shifted to Output Node



4411-0174\_0064

Once the charge within all of the Shift Register's pixels has been shifted out and digitized, the next row of charge is shifted vertically into the shift register where it is then shifted, pixel by pixel into the Output Node and digitized. This process continues until all accumulated charge in the CCD has been shifted out and digitized. The result is zero charge stored within each CCD pixel and the array is immediately ready for the next exposure.

## 9.1 Calculating Image Acquisition/Readout Time

The total time required to acquire and readout a full frame of data at full resolution is calculated as follows:

$$(1) \quad t_{FF} = t_R + t_{exp} + t_c$$

where:

- $t_R$  is the CCD readout time;  
Refer to [Section 9.1.1, CCD Readout Time](#), for additional information.
- $t_{exp}$  is the exposure time;  
This is a user-defined value, and is configured in LightField on the Common Acquisition Setting expander.
- $t_c$  is the sum of the shutter opening and closing compensation times.

### 9.1.1 CCD Readout Time

LightField automatically calculates the readout time, including Region of Interest (ROI) operations and an approximation of overhead times.

CCD readout time is approximately calculated as follows:

$$(2) \quad t_R = (N_y \times t_i) + (N_y \times N_x \times t_{pix}) + overheads$$

where:

- $t_R$  is the readout time;
- $N_y$  is the number of rows within the CCD;
- $N_x$  is the number of columns within the CCD;  
For example, if a CCD is 2048 columns wide,  $N_x = 1024$ .
- $t_i$  is the time needed to shift one line into the shift register;
- $t_{pix}$  is the time needed to process one pixel.

A subsection of the CCD can be read out at full resolution, sometimes dramatically increasing the readout rate while retaining the highest resolution in the Region of Interest (ROI). To approximate the readout rate of an ROI, in Equation 2 substitute the x dimension of the ROI in place of the dimensions of the full CCD. Some overhead time, however, is required to read out and discard the unwanted pixels.

## 9.2 Experiment Timing


For many experiments, the acquisition of quality/useful data is dependent on precise synchronization with external experiment events. Using the EXT SYNC input on the rear of Lansis, externally-generated trigger pulses can be used to control:

- Shutter operation;
- Data readout.

This section describes how to configure Lansis to use incoming trigger pulses to precisely control experiment synchronization in Full Frame mode.

Figure 9-4 illustrates a typical Trigger In expander.

Figure 9-4: Trigger In Expander



4411-0174\_0030

Depending on the specific experiment, two parameters are used to configure Lansis's response to an incoming trigger pulse:

- [Trigger Response](#);
- [Trigger Determined By](#).

Information about the configuration and use of each of these parameters is described in the following sections.

## 9.2.1 Trigger Response

The Trigger Response parameter defines how, upon receipt of an incoming trigger pulse, Lansis reads out data that have been acquired.

Supported Trigger Response modes are:

- No Response;
- Start on Single Trigger;
- Readout Per Trigger;
- Expose During Trigger Pulse.

The following sections describe each of these options and how experiment synchronization is impacted. Within the following sections, the following symbols may be used:

- $t_{\text{exp}}$  = exposure time;
- $t_{\text{O}}$  = shutter opening delay;
- $t_{\text{C}}$  = shutter closing delay;
- $t_{\text{R}}$  = data readout time.

### 9.2.1.1 No Response

When No Response is selected, incoming trigger pulses are ignored. This mode is typically used for experiments incorporating a constant light source (e.g., a CW laser, DC lamp.) Other experiments that can use this mode are high repetition studies where the number of light impulses occurring during a single shutter cycle is so large that the light source appears to be a continuously illuminated source.

Supported Shutter modes are:

- Normal;
- Always Closed;
- Always Open.

The following sections describe how each of these modes impacts experiment timing.


**NOTE:**

In all timing diagrams, Trigger Determined By Rising Edge is illustrated.

### 9.2.1.1.1 Normal

Figure 9-5 illustrates the timing diagram for No Response mode combined with Normal shutter mode.

Figure 9-5: Full Frame Timing Diagram: No Response, Normal




4411-0174\_0031

### 9.2.1.1.2 Always Closed

Figure 9-6 illustrates the timing diagram for No Response mode combined with Always Closed shutter mode.

Figure 9-6: Full Frame Timing Diagram: No Response, Always Closed




4411-0174\_0008

This mode is primarily used when acquiring a dark reference file. However, Lansis typically runs at such a low temperature that a dark reference file is often unnecessary.

### 9.2.1.1.3 Always Open

Figure 9-7 illustrates the timing diagram for No Response mode combined with Always Open shutter mode.

Figure 9-7: Full Frame Timing Diagram: No Response, Always Open



This mode is typically used when an experiment does not support waiting for the shutter to open. Ideally, the only light generated by the experiment is the signal of interest (e.g., a dark chamber with a spark discharge.)

### 9.2.1.2 Start on Single Trigger

Begins the experiment when the trigger is received and the system executes all programmed events.

Supported Shutter modes are:

- Normal;
- Always Closed;
- Always Open;
- Open Before Trigger.

The following sections describe how each of these modes impacts experiment timing.



#### NOTES:

1. When Clean Until Trigger is enabled, an inherent jitter is introduced into the system. Refer to [Section 7.3.1, Clean Until Trigger](#), on page 41 for additional information.
2. In all timing diagrams, Trigger Determined By Rising Edge is illustrated.

#### 9.2.1.2.1 Normal


[Figure 9-8](#) illustrates the timing diagram for Start on Single Trigger mode combined with Normal shutter mode.



#### NOTE:

In [Figure 9-8](#), jitter associated with the cleaning pattern is omitted for clarity. See [Figure 7-3](#) for jitter information.

**Figure 9-8: Full Frame Timing Diagram: Start on Single Trigger, Normal**






### 9.2.1.2.2 Always Closed

Figure 9-9 illustrates the timing diagram for Start on Single Trigger mode combined with Always Closed shutter mode.

Figure 9-9: Full Frame Timing Diagram: Start on Single Trigger, Always Closed




This mode is primarily used when acquiring a dark reference file. However, Lansis typically runs at such a low temperature that a dark reference file is often unnecessary.

### 9.2.1.2.3 Always Open

Figure 9-10 illustrates the timing diagram for Start on Single Trigger mode combined with Always Open shutter mode.

Figure 9-10: Full Frame Timing Diagram: Start on Single Trigger, Always Open



This mode is typically used when an experiment does not support waiting for the shutter to open. Ideally, the only light generated by the experiment is the signal of interest (e.g., a dark chamber with a spark discharge.)

#### 9.2.1.2.4 Open Before Trigger


Figure 9-11 illustrates the timing diagram for Start On Single Trigger mode combined with Open Before Trigger shutter mode.



#### NOTE:

In Figure 9-11, jitter associated with the cleaning pattern is omitted for clarity. See Figure 7-4 for jitter information.

Figure 9-11: Full Frame Timing Diagram: Start on Single Trigger, Open Before Trigger



#### 9.2.1.3 Readout Per Trigger

With Readout Per Trigger, all exposures are synchronized with an incoming trigger pulse. Synchronization occurs on either the Rising Edge or Falling Edge of the trigger pulse which is configured using the Trigger Determined By parameter on the Trigger In expander. Refer to Section 9.2.2, Trigger Determined By, on page 65 for complete information.

Supported Shutter modes are:

- Normal;
- Always Closed;
- Always Open;
- Open Before Trigger.

The effect each of these shutter modes has on experiment synchronization is described in the following sections.




#### NOTES:

1. When Clean Until Trigger is enabled, an inherent jitter is introduced into the system. Refer to Section 7.3.1, Clean Until Trigger, on page 41 for additional information.
2. In all timing diagrams, Trigger Determined By Rising Edge is illustrated.

#### 9.2.1.3.1 Normal

Figure 9-12 illustrates the timing diagram for Readout Per Trigger mode combined with Normal shutter mode.


**Figure 9-12: Full Frame Timing Diagram: Readout Per Trigger, Normal**



### 9.2.1.3.2 Always Closed

Figure 9-13 illustrates the timing diagram for Readout Per Trigger mode combined with Always Closed shutter mode.

Figure 9-13: Full Frame Timing Diagram: Readout Per Trigger, Always Closed




This mode is primarily used when acquiring a dark reference file. However, Lansis typically runs at such a low temperature that the dark reference file is typically unnecessary.

### 9.2.1.3.3 Always Open

Figure 9-14 illustrates the timing diagram for Readout Per Trigger mode combined with Always Open shutter mode.

Figure 9-14: Full Frame Timing Diagram: Readout Per Trigger, Always Open




This mode is typically used when an experiment does not support waiting for the shutter to open. Ideally, the only light generated by the experiment is the signal of interest (e.g., a dark chamber with a spark discharge.) **Open Before Trigger** is closely related.

### 9.2.1.3.4 Open Before Trigger

Figure 9-15 illustrates the timing diagram for Readout Per Trigger mode combined with Open Before Trigger shutter mode.

Figure 9-15: Full Frame Timing Diagram: Readout Per Trigger, Open Before Trigger



When Open Before Trigger is selected, the active shutter is partially synchronized with the experiment.

Upon arrival of the first trigger pulse at the SYNC connector, the shutter:

- Remains open for the configured exposure period;
- Closes;
- The CCD is read out.

Once data readout is complete, the active shutter reopens and waits for the next frame.

Open Before Trigger is most useful when the time between the trigger pulse and the desired incoming signal is not sufficiently long enough to allow the active shutter to completely open.

Unfortunately, this mode exposes the CCD to ambient light while the shutter is open between frames. If the ambient light is constant, triggers occur at regular intervals, and continuous cleaning is enabled, then, depending on the specific cleaning pattern used, most of the average background light can be removed, except for noise attributed to jitter (refer to [Section 7.3.1, Clean Until Trigger](#), on page 41.) If a short cleaning pattern is used, not all ambient light will be removed since only a portion of the CCD will be cleaned on each pattern, and the result is a blurry residual.

Additionally, shot noise of the background cannot be subtracted which may be significant. Consider applications using High gain (i.e.,  $1e^-/ADU$ ), and ambient equals 10% of full scale (i.e.,  $6500 e^-$ .) The shot noise is then approximately 80 ADU, which is significantly greater than the camera's read noise.

As with Normal shutter mode, accurate background subtraction may not be possible for the first frame.

In addition to signal from ambient light, dark charge accumulates during the wait time. Any variation in the external sync frequency also affects the amount of dark charge, even if light is not falling on the CCD during this time.

### 9.2.1.4 Expose During Trigger Pulse

Controls when exposure begins and ends.

Supported Shutter modes are:

- Normal;
- Always Closed;
- Always Open;
- Open Before Trigger.

The effect each of these shutter modes has on experiment synchronization is described in the following sections. In all timing diagrams, Trigger Determined By Rising Edge is illustrated.




#### NOTES:

1. When Clean Until Trigger is enabled, an inherent jitter is introduced into the system Refer to [Section 7.3.1, Clean Until Trigger](#), on page 41 for additional information.
2. In all timing diagrams, Trigger Determined By Rising Edge is illustrated.

#### 9.2.1.4.1 Normal

[Figure 9-16](#) illustrates the timing diagram for Expose During Trigger mode combined with Normal shutter mode.


**Figure 9-16: Full Frame Timing Diagram: Expose During Trigger, Normal**



### 9.2.1.4.2 Always Closed

Figure 9-17 illustrates the timing diagram for Expose During Trigger Pulse mode combined with Always Closed shutter mode.

Figure 9-17: Timing Diagram: Expose During Trigger, Always Closed




This mode is primarily used when acquiring a dark reference file. However, Lansis typically runs at such a low temperature that the dark reference file is typically unnecessary.

### 9.2.1.4.3 Always Open

Figure 9-18 illustrates the timing diagram for Expose During Trigger Pulse mode combined with Always Open shutter mode.

Figure 9-18: Timing Diagram: Expose During Trigger, Always Open




This mode is typically used when an experiment does not support waiting for the shutter to open. Ideally, the only light generated by the experiment is the signal of interest (e.g., a dark chamber with a spark discharge.) [Open Before Trigger](#) is closely related.

#### 9.2.1.4.4 Open Before Trigger

Figure 9-19 illustrates the timing diagram for Expose During Trigger Pulse mode combined with Open Before Trigger shutter mode.

Figure 9-19: Timing Diagram: Expose During Trigger, Open Before Trigger





## 9.2.2 Trigger Determined By

When using an external trigger to initiate a readout, Lansis can be configured to respond to:

- The rising edge of the incoming trigger pulse;
- The falling edge of the incoming trigger pulse.

The Trigger Determined By parameter configures this behavior. Valid values are:

- Rising Edge;

Lansis responds to the rising edge of incoming trigger pulses. Depending on the specific system configuration, one or more subsequent trigger pulses may be required (or ignored) by the system. Refer to [Section 9.2.1, Trigger Response](#), on page 53 for complete information about configuring incoming trigger responses.

When using [Expose During Trigger Pulse](#) mode:

- Exposure begins on the Rising Edge of each incoming trigger pulse;
- Exposure ends on the Falling Edge of the respective trigger pulse.

The timing diagrams in [Section 9.2.1.4, Expose During Trigger Pulse](#), on page 62 illustrate this.



---

**NOTE:**

If Lansis is busy when a subsequent trigger pulse is received, the trigger is ignored.

---

- Falling Edge

Lansis responds to the falling edge of incoming trigger pulses. Depending on the specific system configuration, one or more subsequent trigger pulses may be required (or ignored) by the system. Refer to [Section 9.2.1, Trigger Response](#), on page 53 for complete information about configuring incoming trigger responses.

When using [Expose During Trigger Pulse](#) mode:

- Exposure begins on the Falling Edge of each incoming trigger pulse;
- Exposure ends on the Rising Edge of the respective trigger pulse.



---

**NOTE:**


If Lansis is busy when a subsequent trigger pulse is received, the trigger is ignored.

---

## 9.3 Trigger Out

In addition to being able to synchronize Lansis with an experiment, additional equipment can be synchronized using the Trigger Out connector on the rear of Lansis. Two trigger out pulses are configured on the Trigger Out expander, shown in [Figure 9-20](#).

Figure 9-20: Typical Trigger Out Expander



4411-0174\_0034

The following options are available for each of the two output signals:

- Acquiring;  
The associated output signal is high when Lansis is acquiring or ready to receive the first trigger.
- Always High;  
The associated output signal is always high.
- Exposing;  
The associated output signal is high when the sensor is exposed as configured within LightField.
- Reading Out;  
The associated output signal when data are being read out of the sensor.
- Shutter Open;  
The associated output signal is high when the shutter is open.
- Waiting for Trigger.  
The associated output signal when Lansis is waiting for an incoming trigger.

Each of these options can also be inverted to create active low signals using the Invert Output Signal option.

# Chapter 10: Binning

---

Binning is the process of summing data from adjacent pixels to form a single pixel, often called a Super Pixel. Binning can be accomplished in one of two ways:

- Hardware;
- Software.

Rectangular groups of pixels of any size may be binned together subject to some hardware and software limitations.


## 10.1 Hardware Binning

Hardware binning is performed on the CCD array before the signal is read out of the output amplifier. For signal levels that are readout noise limited this method improves S/N ratio linearly with the number of pixels grouped together. For signals large enough to render the camera photon shot noise limited, the S/N ratio improvement is roughly proportional to the square-root of the number of pixels binned.

Binning reduces readout time and the burden on computer memory, but at the expense of resolution. Since shift register pixels typically hold only twice as much charge as image pixels, the binning of large sections may result in saturation and “blooming”, or spilling of charge back into the image area.

[Figure 10-1](#) illustrates an example of 2 x 2 binning. Each pixel of the image displayed by the software represents 4 pixels of the CCD array. Rectangular bins of any size are possible.

Figure 10-1: 2 × 2 Binning



A3  
A4  
A5  
A6

B3  
B4  
B5  
B6

C3  
C4  
C5  
C6

D3  
D4  
D5  
D6

A3  
A4  
A5  
A6

B3  
B4  
B5  
B6

C3  
C4  
C5  
C6

D3  
D4  
D5  
D6

A3  
A4  
A5  
A6

B3  
B4  
B5  
B6

C3  
C4  
C5  
C6

D3  
D4  
D5  
D6

4411-0174\_0040

## 10.2 Software Binning

One limitation of hardware binning is that the shift register pixels and the output node are typically only 2-3 times the size of imaging pixels. Consequently, if the total charge binned together exceeds the capacity of the shift register or output node, the data will be corrupted.

This restriction strongly limits the number of pixels that may be binned in cases where there is a small signal superimposed on a large background, such as signals with a large fluorescence. Ideally, one would like to bin many pixels to increase the S/N ratio of the weak peaks but this cannot be done because the fluorescence would quickly saturate the CCD.

The solution is to perform the binning in software. Limited hardware binning may be used when reading out the CCD. Additional binning is accomplished in software, producing a result that represents many more photons than was possible using hardware binning.

Software averaging can improve the S/N ratio by as much as the square-root of the number of pixels binned. Unfortunately, with a high number of pixels binned, i.e., above 100, camera 1/f noise may reduce the actual S/N ratio to slightly below this theoretical value. Also, if the light source used is photon-flicker limited rather than photon shot-noise limited, this theoretical signal improvement cannot be fully realized. Again, background subtraction from the raw data is necessary.

This technique is also useful in high light level experiments (e.g., absorbance spectroscopy,) where the camera is again photon shot-noise limited. Summing multiple pixels in software corresponds to collecting more photons, and results in a better S/N ratio in the measurement.

*This page is intentionally blank.*


# Chapter 11: Shutter Control

This chapter provides information about the configuration and control of an external shutter that is connected to a Lansis camera. Lansis cameras are not supplied with a driver for internal or external shutter, and an external shutter driver is required to operate the shutter. The timing information below will still apply, but the software then controls the timing of the shutter pulse applied to a logic out MCX connector, which should in turn be connected to the external shutter driver. You should select the "Shutter Open" signal in the Trigger Out expander.

## 11.1 Configuration

Shutter information is configured within LightField on the Shutter expander. [Figure 11-1](#) illustrates a typical Shutter expander.

Figure 11-1: Typical Shutter Expander:



### 11.1.1 Mode

This parameter determines the shutter's opening and closing behavior during an experiment.

Depending upon the specific Trigger In ► Trigger Response that has been selected, supported Modes are:

- Normal;  
The shutter opens for exposure, and closes when complete.
- Always Closed;  
When selected, the shutter closes and remains in the closed position.  
This mode is primarily used when acquiring a dark reference file. However, Lansis typically runs at such a low temperature that a dark reference file is often unnecessary.
- Always Open;  
When selected, the shutter opens and remains open for the experiment duration.  
This mode is typically used when an experiment does not support waiting for the shutter to open. Ideally, the only light generated by the experiment is the signal of interest (e.g., a dark chamber with a spark discharge.)

- **Open Before Trigger**  
When selected, the shutter opens as soon as Lansis is ready to receive an External Sync pulse. This is required if the time between the External Sync pulse and the event is less than the time it takes the shutter to open which is a typically a few milliseconds.

### **11.1.2 Opening Delay**

Specifies, in milliseconds, the length of time Lansis is to wait for the shutter to open.

Valid values are [0 ... 1000] ms, inclusive, in 1 ms increments.

### **11.1.3 Closing Delay**

Specifies, in milliseconds, the length of time Lansis is to wait for the shutter to close.

Valid values are [0 ... 1000] ms, inclusive, in 1 ms increments.



## 11.2 Using an External Shutter



### **WARNING! RISK OF ELECTRIC SHOCK!**

---

When using a shutter not supplied by Teledyne Princeton Instruments, it is the responsibility of each user to assure all connections, other than ground, to the external shutter are insulated to prevent accidental contact by personnel.

Although not mandatory, Teledyne Princeton Instruments strongly recommends the use of shielded cable with braided shield between the shutter and the shutter connector. The shield braid should be connected to the shell of the Hirose connector in accordance with Hirose's recommendations for connecting shield grounds.

Usually, it is best to connect the shield ground to the frame of the external shutter, but this may vary with the circumstances. Although the use of a braided shield typically reduces the susceptibility of a system to electromagnetic disturbances, this protection is not guaranteed in all applications.

---



### **CAUTION!**

---

In order to prevent potential permanent damage to either the Lansis camera and/or the shutter, always contact Teledyne Princeton Instruments before connecting an external shutter. Refer to [Contact Information](#) on page 100 for complete information.

---

*This page is intentionally blank.*

# Appendix A: Technical Specifications



## CAUTION!

All specifications are subject to change.

This appendix provides some technical information and specifications for Lansis cameras and optional accessories. Additional information may be found on data sheets available on the Teledyne Princeton Instruments website ([www.princetoninstruments.com](http://www.princetoninstruments.com)).

## A.1 System Dimensions and Weight

Refer to [Table A-1](#) for system dimensions and weight.

**Table A-1: General Camera Specifications**

Dimension	Specification
Length	5.913 in [15.019 cm]
Width	4.863 in [12.35 cm]
Height	4.863 in [12.35 cm]
Weight	6.83 lbs [3100 g]

**Table A-2: General Power Supply**

Dimension	Specification
Length	1.38 in [35 mm]
Width	2.64 in [67 mm]
Height	6.57 in [167 mm]

### A.1.1 Vacuum Window

SI-UV fused-silica quartz (0.125"/3.17 mm thick)

## A.2 Camera Specifications

Refer to [Table A-3](#) for CCD specifications for Lansis detectors.

**Table A-3: Lansis CCD Specifications<sup>a</sup> (Sheet 1 of 2)**

Specification	A-261BR	A-261BRX	W-261BR	W-261BRX
CCD	TE2V CCD261			
Resolution	2048 x 263			
Pixel Size	15 $\mu\text{m}$ x 15 $\mu\text{m}$			

**Table A-3: Lansis CCD Specifications<sup>a</sup> (Sheet 2 of 2)**

Specification	A-261BR	A-261BRX	W-261BR	W-261BRX
Imaging Area	30.72 mm x 3.96 mm			
Spectra/s <sup>b</sup> (full vertical binning)	~236 fps			
ADC Speed/16 bits	100 kHz, 2 MHz			
Vertical Shift Rate <sup>c</sup>	7.41 $\mu$ sec/row			
Read Noise	5 e <sup>-</sup> rms @ 100 kHz (typical)			
Optical Focal Distance	TBD			

- a. Specifications are valid as of the publication date of this manual. For up-to-date specifications, refer to the Lansis data sheets available for download from [www.princetoninstruments.com](http://www.princetoninstruments.com).
- b. @ 2 MHz full vertical binning
- c. Software programmable

### A.2.1 Thermal Characteristics

Refer to [Table A-4](#) for specific thermal information.

**Table A-4: Default Operating Temperature**

Specification	Cooling Medium	Specification
Default Operating Temperature	Air	-60°C
	Liquid <sup>a</sup>	-65 °C
Deepest Cooling Temperature	Air	-65 °C
	Liquid <sup>a</sup>	-70 °C
Precision	–	±0.05 °C

- a. External chiller required.

## A.3 Power Specifications

All voltages required by Lansis cameras are generated and delivered by an external power supply included with each Lansis camera using the supplied cables.

Power to the camera is switched on and off using the toggle switch on the external power supply.



### **WARNING!**

In case of a fire or other emergency, immediately remove the power supply's AC plug from the wall receptacle.



### **CAUTION!**

Use of a power supply other than that provided with the Lansis camera will void the camera warranty. For specific power supply requirements, contact Teledyne Princeton Instruments. Refer to [Contact Information](#) on page 100 for complete information.

Refer to [Table A-5](#) for power specifications for the external Lansis power supplies.

**Table A-5: Power Specifications**

Parameter	Specification	Units
Input Voltage (nominal)	100 – 240	V <sub>AC</sub>
Input Frequency (nominal)	47 – 63	Hz
Input Power (maximum)	160	W

## A.4 Environmental Specifications

Refer to [Table A-6](#) for environmental specifications.

**Table A-6: Lansis Environmental Specifications**

Parameter	Specification		
	Minimum	Nominal	Maximum
Storage Temperature	-20°C		+55°C
Operating Temperature	+5°C		+30°C
Operating Ambient Relative Humidity	<80% (non-condensing)		
Operating Ambient Temperature	0°C		+25°C



### NOTE:

Cooling performance may degrade if the room temperature is above +23°C.

### A.4.1 Ventilation

A minimum of 1 inch (2.54 cm) clearance is required around all vents on the Lansis camera.

Where Lansis is within an enclosure, >30 cfm air circulation and heat dissipation of 120 W is required.

## A.5 External Chiller Specifications



### CAUTION!

Never set the coolant temperature below the dew point.

Refer to [Table A-7](#) for external chiller specifications.

**Table A-7: External Chiller Specifications**

Parameter	Specification
Flow	0.8 gal/minute minimum
Pressure	30.0 psi maximum
Coolant	50/50 Glycol Water
Coolant Temperature	20°C non condensing

## A.6 Minimum Host Computer Specifications



### NOTE:

Computers and operating systems experience frequent updates. Therefore, the following sections are intended to provide minimum system requirements for operating a Lansis camera.

A faster computer with 5 GB or larger memory (RAM) will greatly enhance the software performance during live mode operations.

Contact the factory to determine specific requirements.

---

The minimum system requirements for LightField are:

- Windows 7/8/10 (64-bit)
- 2 GHz dual core processor
- 4 GB RAM (or greater)
- CD-ROM drive
- Super VGA monitor and graphics card supporting at least 65535 colors with at least 128 MB of memory. Memory requirement is dependent on desired display resolution.
- Hard disk with a minimum of 1 GB available for installation. Additional space is required for data storage: the amount of space required depends on the number and size of images/spectra collected. Disk level compression programs are not recommended. Drive speed of 10,000 RPM recommended.
- Mouse or other pointing device.


*This page is intentionally blank.*



# Appendix B: Outline Drawings


This appendix provides outline drawings for the Lansis Camera System.

Figure B-1: Outline Drawing: Lansis Camera Mount



4411-0174\_0107

**Figure B-2: Outline Drawing: Lansis Spectrometer Mount**



4411-0174\_0107


# Appendix C: Drain Coolant from Lansis

---

This appendix provides information necessary to safely drain coolant from within the Lansis camera body.

1. Place Lansis camera body on a flat, secure surface with the two coolant fittings facing down and positioned over a container that will collect the coolant as it is draining. See [Figure C-1](#).


**Figure C-1: Positioning Lansis to Drain Coolant**



4411-0174\_0128

2. Using two small screwdrivers, Allen wrenches, or similar tools, carefully depress the center piston on each coolant fitting simultaneously. See [Figure C-2](#)

**Figure C-2: Depressing Coolant Fittings' Center Pistons**



**Carefully depress using two small screwdrivers, allen wrenches, or similar tool.**

This will allow the coolant to drain out of the Lansis and into the collection bucket.



**NOTE:**

It is not necessary to drain every drop of coolant out of the Lansis, but always drain as much as reasonably possible.

3. When the coolant has drained, release the pistons and verify that they return to their normal, closed position.
4. Use a soft cloth to wipe off any stray coolant from the Lansis and/or surrounding areas.
5. If available, install rubber fitting covers onto the coolant fittings prior to shipping or storing the Lansis.
6. Dispose of drained coolant according to local standards/requirements.

EAR99 Technology Subject to Restrictions Contained on the Cover Page.

# Appendix D: Custom Modes

Custom Sensor and Custom Timing are standard with LightField, although both are sensor and readout mode-dependent. These modes enable data acquisition at the fastest possible rates for the camera. Custom Sensor allows the apparent size of the CCD array to be reduced, while Custom Timing allows a faster vertical shift time to be selected.

## D.1 Custom Sensor

In addition to Binning and ROI (previously discussed in the manual), the Custom Sensor feature can be used to reduce Readout Time.




### NOTE:

Teledyne Princeton Instruments does not encourage users to change these parameter settings. For most applications, the default settings will give the best results. There are limitations to this technique, such as requiring the careful positioning of the light on the active section of the CCD and darkness on the rest of it.

Custom Sensor redefines the size of the CCD's active area via software. Unlike setting a smaller Region of Interest (ROI), which also involves reading out fewer pixels, this mode does not incur overhead from discarding or skipping the rest of the rows. And, unlike both Binning and ROI, Custom Sensor relies on a form of array masking to ensure that no light falls outside the currently set active area.

The Custom Sensor pane, illustrated in [Figure D-1](#), is accessed by opening the Sensor expander and clicking on the Custom Sensor button.

Figure D-1: Custom Sensor Pane



4411-0174\_0043

By changing the values in the Active fields, the image acquisition speed can be increased by reducing the size of the active area in the definition. The result will be faster, but lower resolution, data acquisition. Operating in this mode would ordinarily require that the chip be masked so that only the reduced active area is exposed. This will prevent unwanted charge from spilling into the active area or being transferred to the shift register.

With a high performance imaging spectrometer, such as Teledyne Princeton Instruments' IsoPlane family, this can be done by masking at the input of the spectrometer. This has the advantage of preventing the need to add masking inside the camera. Older imaging spectrometers can be used this way, but the results will not be as good as with an IsoPlane.

By default, if there are no Pre-Dummy rows, the serial register will be cleared before rows are shifted.

**NOTE:**

The Clean Serial Register function only appears in the Sensor Cleaning pane when the Inactive Area Top Margin is 0 rows and Pixel Bias Correction is off. De-select the check box to deactivate the serial register cleaning.

### D.1.1 Custom Timing


**NOTE:**

Custom Timing is standard within LightField for full frame CCD cameras.

Custom Timing is configured using the Custom Sensor button on the Sensor expander. Vertical shift rate is the time required to shift one row into the serial register. The smaller the value, the faster charge will be shifted up one row at a time toward the serial register.

In the Custom Timing panel, vertical shift rate is configured by selecting the desired rate from the pull-down menu, as illustrated in [Figure D-2](#).

**Figure D-2: Custom Timing**



4411-0174\_0008

# Appendix E: Troubleshooting



## **WARNING!**

Do not attach or remove any cables while the camera system is powered on.

**Table E-1: Troubleshooting Index by Error/Fault Description**

Error/Fault	Refer to...
Power Supply Switch in On Position, But Power LED Extinguished	page 89
Overexposed CCD	page 90
Camera Stops Working	page 90
Temperature Lock Cannot be Achieved or Maintained	page 91
Camera Loses Temperature Lock	page 92
Gradual Deterioration of Cooling Capability	page 92
External Coolant Circulator: Low Coolant (Air in the Hoses)	page 92
LightField Faults/Errors	page 93
Devices Missing	page 93
Device is Occupied	page 94
Acquisition Started but Viewer Display Does Not Update	page 94

## **E.1 General Camera Faults/Errors**

This section provides information about troubleshooting general camera faults and errors.

### **E.1.1 Connection Failure or Logic Power Supply Overcurrent**

This failure is indicated by the red FAULT LED on the rear of the Lansis power supply being ON (i.e., does not flash.)

Examine the primary power cable and repair/replace as necessary.

Verify all connectors are fully seated at both ends and the jack screws are secured.

### **E.1.2 Power Supply Switch in On Position, But Power LED Extinguished**

If the power supply has been turned on (i.e., power switch is in the I position,) but the green POWER LED does not illuminate indicating the power supply is operational:

- Turn the power switch off (i.e., the 0 position);
- Wait approximately five seconds;
- Turn the power switch on (i.e., the I position)

If the green POWER LED still remains extinguished, the power supply must be returned to Teledyne Princeton Instruments for repair. Refer to [Contact Information](#) on page 100 for complete information.

### E.1.3 Overexposed CCD

It takes an enormous power density to damage the CCD at room temperature or cooled. It can be done (e.g., imaging the sun with a large aperture lens, so that the silicon is heated.)

What can happen if the CCD is cold is a temporary elevation of the dark current. If the CCD is over-saturated, electrons pool in the substrate and are clocked out slowly (over hours or even days). To clear the charge, it is best to return the CCD to room temperature, then run the CCD (set short exposure, then click run in LF) while it is cooled down to normal operating temperature.

### E.1.4 Camera Stops Working

Problems with the host computer system or software may have side effects that appear to be hardware problems. If you are sure the problem is in the camera system hardware, perform these preliminary system checks:

- Examine the STATUS LED on the rear panel of the camera. See [Table 2-2](#) for decoding STATUS LED states.
- Turn off all AC power.
- Verify that all cables are securely fastened.
- Turn the system back on.

If the system still does not respond, contact Customer Support. Refer to [Contact Information](#) on page 100 for complete information.



## E.2 Cooling Faults/Errors

This section provides recommended troubleshooting guidelines for cooling-related issues.

### E.2.1 Temperature Lock Cannot be Achieved or Maintained



#### CAUTION!

The most probable cause of a failure to lock is the setpoint has been programmed for a temperature lower than Lansis can achieve. Return the temperature setpoint to the default value. Protracted operation with an unachievable setpoint may accelerate TEC failure.

Possible causes for not being able to achieve or maintain lock include:

- Ambient temperature greater than +25°C. If the ambient temperature is greater than +25°C, either cool the camera's environment or raise the set temperature.
- Airflow through the camera and/or circulator is obstructed. The camera needs to have approximately 2" [50.8 mm] of clearance around the vented covers. If there is an enclosure involved, the enclosure needs to have unrestricted flow to an open environment. The camera vents its heat out the vents near the nose. The air intake is near the rear of the camera.
- If using liquid cooling, additional causes for a failure to achieve and/or maintain lock may include:
  - A hose is kinked. Unkink the hose.
  - Coolant level is low. Add coolant. Refer to [Section E.2.4, External Coolant Circulator: Low Coolant \(Air in the Hoses\)](#), on page 92.
  - There may be air in the hoses. Add coolant. Refer to [Section E.2.4, External Coolant Circulator: Low Coolant \(Air in the Hoses\)](#), on page 92.
  - Circulator pump is not working. If you do not hear the pump running when the external chiller is powered on, turn off the circulator and contact Customer Support. Refer to [Contact Information](#) on page 100 for complete information.
  - The circulator is higher than the camera. Reposition the circulator so that it is 6" [150 mm] or more below the camera. The vertical distance should not exceed 10 feet [3 m]. Typically, the camera is at table height and the circulator is on the floor.
- The camera's internal temperature may be too high which may occur if the operating environment is particularly warm or when attempting to operate at a temperature colder than the specified limit. Lansis cameras are equipped with a thermal-protection switch that shuts the cooler circuits down if the internal temperature exceeds a preset limit. Typically, camera operation is restored automatically after approximately ten minutes. Although the thermo-protection switch will protect the camera, it is nevertheless advised to power down the camera and correct the operating conditions that caused the thermal-overload to occur.
- The camera vacuum has deteriorated and needs to be refreshed. Contact Customer Support. Refer to [Contact Information](#) on page 100 for complete information.

### **E.2.2 Camera Loses Temperature Lock**

The internal temperature of the camera is too high. This might occur when the operating environment is particularly warm or when attempting to operate at a temperature colder than the specified limit. If this happens, an internal thermal overload switch will disable the cooler circuits to protect them. Typically, camera operation is restored in about ten minutes. Although the thermal overload switch will protect the camera, users are advised to power down and correct the operating conditions that caused the thermal overload to occur.

Additionally, repeated cycling can reduce the lifetime of the thermoelectric cooling system.

### **E.2.3 Gradual Deterioration of Cooling Capability**

While unlikely with the Lansis camera (guaranteed permanent vacuum for the life of the camera), if a gradual deterioration of the cooling capability is observed, there may be a gradual deterioration of the camera's vacuum. This can affect temperature performance such that it may be impossible to achieve temperature lock at the lowest temperatures. In the kind of applications for which cooled CCD cameras are so well suited, it is highly desirable to maintain the system's lowest temperature performance because lower temperatures result in lower thermal noise and better the signal-to-noise ratio. Contact the factory to make arrangements for returning the camera to the support facility. Refer to [Contact Information](#) on page 100 for complete information.

### **E.2.4 External Chiller: Low Coolant (Air in the Hoses)**

Most circulators and/or chillers will not operate properly if there is insufficient coolant in their tank or if excess air gets into the system. It is advised that you check the manual for the circulator/chiller in use.


## E.3 LightField Faults/Errors

This section provides information about troubleshooting problems that may occur with LightField.

### E.3.1 Devices Missing

When LightField is started, it looks for devices that are powered on and connected via a communications interface to the host computer. If it cannot find a device that was used in the last experiment, the dialog shown in [Figure E-1](#) is displayed while LightField continues searching for the missing device.

Figure E-1: Devices Missing Dialog



4411-0174\_0068


Perform the following steps to try to resolve this fault:

- Verify the device is connected and powered on.  
If the device is connected but turned off, switch it on.  
LightField should now find the device. If it does not, cancel the load and restart LightField.
- Cancel the loading of the experiment.  
Canceling an experiment's loading means that the last used experiment will not be loaded automatically when LightField opens. However, the experiment may be loaded after all devices are available, a new experiment design can be started, or a different experiment can be loaded which uses the devices that are available.

### E.3.2 Device is Occupied

Although multiple instances of LightField can be running at the same time, any device that is currently in use by one instance of LightField will be shown within the Available Devices area of all other instances as Occupied. See [Figure E-2](#).

Figure E-2: Typical LightField Occupied Device




4411-0174\_0071

To have a device become available to the current instance of LightField, either remove it from the Experiment Devices area in the other instance or close the instance that is using the device.

### E.3.3 Acquisition Started but Viewer Display Does Not Update


Live data being acquired in either Preview or Acquire mode are displayed in a Data Viewer tab on the Experiment workspace. If the active Data Viewer's display is not being updated and data acquisition is occurring, determine if there is a filename displayed within the active Data View tab. See [Figure E-3](#).

Figure E-3: Typical Acquisition Display




4411-0174\_0069

When a filename is listed, it indicates that the data being displayed are static (i.e., data from the indicated file,) and not live data that are currently being acquired.

To return to a live data view, click on the  to the right of the filename to view the pull-down menu and select Live Data. See [Figure E-4](#).

**Figure E-4: Data Viewer Menu**





# Warranty and Service

---

## Limited Warranty

Teledyne Princeton Instruments ("us," "we," "our,") makes the following limited warranties. These limited warranties extend to the original purchaser ("You," "you,") only and no other purchaser or transferee. We have complete control over all warranties and may alter or terminate any or all warranties at any time we deem necessary.

### ***Basic Limited One (1) Year Warranty***

Teledyne Princeton Instruments warrants this product against substantial defects in materials and/or workmanship for a period of up to one (1) year after shipment. During this period, Teledyne Princeton Instruments will repair the product or, at its sole option, repair or replace any defective part without charge to you. You must deliver the entire product to the Teledyne Princeton Instruments factory or, at our option, to a factory-authorized service center. You are responsible for the shipping costs to return the product. International customers should contact their local Teledyne Princeton Instruments authorized representative/distributor for repair information and assistance, or visit our technical support page at [www.princetoninstruments.com](http://www.princetoninstruments.com).

### ***Limited One (1) Year Warranty on Refurbished or Discontinued Products***

Teledyne Princeton Instruments warrants, with the exception of the CCD imaging device (which carries NO WARRANTIES EXPRESS OR IMPLIED,) this product against defects in materials or workmanship for a period of up to one (1) year after shipment. During this period, Teledyne Princeton Instruments will repair or replace, at its sole option, any defective parts, without charge to you. You must deliver the entire product to the Teledyne Princeton Instruments factory or, at our option, a factory-authorized service center. You are responsible for the shipping costs to return the product to Teledyne Princeton Instruments. International customers should contact their local Teledyne Princeton Instruments representative/distributor for repair information and assistance or visit our technical support page at [www.princetoninstruments.com](http://www.princetoninstruments.com).

### ***XP Vacuum Chamber Limited Lifetime Warranty***

Teledyne Princeton Instruments warrants that the cooling performance of the system will meet our specifications over the lifetime of an XP style detector (has all metal seals) or Teledyne Princeton Instruments will, at its sole option, repair or replace any vacuum chamber components necessary to restore the cooling performance back to the original specifications at no cost to the original purchaser. *Any failure to "cool to spec" beyond our Basic (1) year limited warranty from date of shipment, due to a non-vacuum-related component failure (e.g., any components that are electrical/electronic) is NOT covered and carries NO WARRANTIES EXPRESSED OR IMPLIED.* Responsibility for shipping charges is as described above under our Basic Limited One (1) Year Warranty.

### ***Sealed Chamber Integrity Limited 12 Month Warranty***

Teledyne Princeton Instruments warrants the sealed chamber integrity of all our products for a period of twelve (12) months after shipment. If, at anytime within twelve (12) months from the date of delivery, the detector should experience a sealed chamber failure, all parts and labor needed to restore the chamber seal will be covered by us. *Open chamber products carry NO WARRANTY TO THE CCD IMAGING DEVICE, EXPRESSED OR IMPLIED.* Responsibility for shipping charges is as described above under our Basic Limited One (1) Year Warranty.

### ***Vacuum Integrity Limited 12 Month Warranty***

Teledyne Princeton Instruments warrants the vacuum integrity of “Non-XP” style detectors (do not have all metal seals) for a period of up to twelve (12) months from the date of shipment. We warrant that the detector head will maintain the factory-set operating temperature without the requirement for customer pumping. Should the detector experience a Vacuum Integrity failure at anytime within twelve (12) months from the date of delivery all parts and labor needed to restore the vacuum integrity will be covered by us. Responsibility for shipping charges is as described above under our Basic Limited One (1) Year Warranty.

### ***Image Intensifier Detector Limited One Year Warranty***

All image intensifier products are inherently susceptible to Phosphor and/or Photocathode burn (physical damage) when exposed to high intensity light. Teledyne Princeton Instruments warrants, with the exception of image intensifier products that are found to have Phosphor and/or Photocathode burn damage (which carry NO WARRANTIES EXPRESSED OR IMPLIED,) all image intensifier products for a period of one (1) year after shipment. *Refer to additional Limited One (1) year Warranty terms and conditions above, which apply to this warranty.* Responsibility for shipping charges is as described above under our Basic Limited One (1) Year Warranty.

### ***X-Ray Detector Limited One Year Warranty***

Teledyne Princeton Instruments warrants, with the exception of CCD imaging device and fiber optic assembly damage due to X-rays (which carry NO WARRANTIES EXPRESSED OR IMPLIED,) all X-ray products for one (1) year after shipment. *Refer to additional Basic Limited One (1) year Warranty terms and conditions above, which apply to this warranty.* Responsibility for shipping charges is as described above under our Basic Limited One (1) Year Warranty.

### ***Software Limited Warranty***

Teledyne Princeton Instruments warrants all of our manufactured software discs to be free from substantial defects in materials and/or workmanship under normal use for a period of one (1) year from shipment. Teledyne Princeton Instruments does not warrant that the function of the software will meet your requirements or that operation will be uninterrupted or error free. You assume responsibility for selecting the software to achieve your intended results and for the use and results obtained from the software. In addition, during the one (1) year limited warranty. The original purchaser is entitled to receive free version upgrades. Version upgrades supplied free of charge will be in the form of a download from the Internet. Those customers who do not have access to the Internet may obtain the version upgrades on a CDROM from our factory for an incidental shipping and handling charge. *Refer to Item 12 in [Your Responsibility](#) of this warranty for more information.*



## ***Owner's Manual and Troubleshooting***

You should read the owner's manual thoroughly before operating this product. In the unlikely event that you should encounter difficulty operating this product, the owner's manual should be consulted before contacting the Teledyne Princeton Instruments technical support staff or authorized service representative for assistance. If you have consulted the owner's manual and the problem still persists, please contact the Teledyne Princeton Instruments technical support staff or our authorized service representative. *Refer to Item 12 in [Your Responsibility](#) of this warranty for more information.*

## ***Your Responsibility***

The above Limited Warranties are subject to the following terms and conditions:

1. You must retain your bill of sale (invoice) and present it upon request for service and repairs or provide other proof of purchase satisfactory to Teledyne Princeton Instruments.
2. You must notify the Teledyne Princeton Instruments factory service center within (30) days after you have taken delivery of a product or part that you believe to be defective. With the exception of customers who claim a "technical issue" with the operation of the product or part, all invoices must be paid in full in accordance with the terms of sale. Failure to pay invoices when due may result in the interruption and/or cancellation of your one (1) year limited warranty and/or any other warranty, expressed or implied.
3. All warranty service must be made by the Teledyne Princeton Instruments factory or, at our option, an authorized service center.
4. Before products or parts can be returned for service you must contact the Teledyne Princeton Instruments factory and receive a return authorization number (RMA.) Products or parts returned for service without a return authorization evidenced by an RMA will be sent back freight collect.
5. These warranties are effective only if purchased from the Teledyne Princeton Instruments factory or one of our authorized manufacturer's representatives or distributors.
6. Unless specified in the original purchase agreement, Teledyne Princeton Instruments is not responsible for installation, setup, or disassembly at the customer's location.
7. Warranties extend only to defects in materials or workmanship as limited above and do not extend to any product or part which:
  - has been lost or discarded by you;
  - has been damaged as a result of misuse, improper installation, faulty or inadequate maintenance, or failure to follow instructions furnished by us;
  - has had serial numbers removed, altered, defaced, or rendered illegible;
  - has been subjected to improper or unauthorized repair;
  - has been damaged due to fire, flood, radiation, or other "acts of God," or other contingencies beyond the control of Teledyne Princeton Instruments; or
  - is a shutter which is a normal wear item and as such carries a onetime only replacement due to a failure within the original 1 year Manufacturer warranty.
8. After the warranty period has expired, you may contact the Teledyne Princeton Instruments factory or a Teledyne Princeton Instruments-authorized representative for repair information and/or extended warranty plans.
9. Physically damaged units or units that have been modified are not acceptable for repair in or out of warranty and will be returned as received.

10. All warranties implied by state law or non-U.S. laws, including the implied warranties of merchantability and fitness for a particular purpose, are expressly limited to the duration of the limited warranties set forth above. With the exception of any warranties implied by state law or non-U.S. laws, as hereby limited, the forgoing warranty is exclusive and in lieu of all other warranties, guarantees, agreements, and similar obligations of manufacturer or seller with respect to the repair or replacement of any parts. In no event shall Teledyne Princeton Instruments' liability exceed the cost of the repair or replacement of the defective product or part.
11. This limited warranty gives you specific legal rights and you may also have other rights that may vary from state to state and from country to country. Some states and countries do not allow limitations on how long an implied warranty lasts, when an action may be brought, or the exclusion or limitation of incidental or consequential damages, so the above provisions may not apply to you.
12. When contacting us for technical support or service assistance, please refer to the Teledyne Princeton Instruments factory of purchase, contact your authorized Teledyne Princeton Instruments representative or reseller, or visit our technical support page at [www.princetoninstruments.com](http://www.princetoninstruments.com).

## Contact Information

Teledyne Princeton Instruments' manufacturing facility for this product is located at the following address:

Teledyne Princeton Instruments  
3660 Quakerbridge Road  
Trenton, NJ 08619 (USA)

Tel: 1-800-874-9789 / 1-609-587-9797

Customer Support E-mail: [pi.techsupport@teledyne.com](mailto:pi.techsupport@teledyne.com)

Refer to <http://www.princetoninstruments.com/support> for complete support and contact information, including:

- Up-to-date addresses and telephone numbers;
- Software downloads;
- Product manuals;
- Support topics for Teledyne Princeton Instruments' product lines.

All rights reserved. No part of this publication may be reproduced by any means without the written permission of Teledyne Princeton Instruments.

Printed in the United States of America.

Lansis, IntelliCal, IsoPlane, LightField, and PVCAM are registered trademarks of Teledyne Digital Imaging US, Inc..

LabVIEW is a registered trademark of National Instruments, Inc.

LEMO is a registered trademark of INTERLEMO HOLDING SA

Windows is a registered trademark of Microsoft Corporation in the United States and/or other countries.

The information in this publication is believed to be accurate as of the publication release date. However, Teledyne Princeton Instruments does not assume any responsibility for any consequences including any damages resulting from the use thereof. The information contained herein is subject to change without notice. Revision of this publication may be issued to incorporate such change.

*This page is intentionally blank.*



[www.princetoninstruments.com](http://www.princetoninstruments.com)



**TELEDYNE PRINCETON INSTRUMENTS**  
Everywhereyoulook™

Part of the Teledyne Imaging Group

[pi.info@teledyne.com](mailto:pi.info@teledyne.com)

**USA +1 877-474-2286 | Please visit [www.princetoninstruments.com](http://www.princetoninstruments.com) to find the contact information for our global offices.**

**EAR99 Technology Subject to Restrictions Contained on the Cover Page.**